

One Coat **Bond SL**TM

All the resistance you need, in
just one drop!


One Coat Bond SL

One Coat Bond SL is a Total Etch adhesive for dentin and enamel with water based solvent. Its hydrophilic chemical reaction and the balanced composition allow complete penetration into the dentinal tubules, with perfect humidification of the dental structure, thereby promoting homogeneity in the resin layer.

- Compatible with water - Excellent performance on dry and wet dentin
- No loss of product into the environment
- Greater coverage power: just one drop and one layer for perfect adhesion
- Absence of volatile products that could be trapped under the veneer, jeopardizing adhesion
- Does not cause postoperative sensitivity

Furthermore, One Coat Bond SL has the best clinical performance proven in studies conducted by USP Bauru, USP Araraquara, USP Ribeirão Preto, the Federal University of Santa Catarina and Ponta Grossa State University:

Influence of thermocycling on the shear strength of dentin and enamel adhesive systems (500 cycles between 4 °C and 55 °C)



PATRICIA ALEIXO DOS SANTOS, Department of Operative Dentistry, USP/Bauru
PATRICIA PETROMILLI NORDI SASSO GARCIA, Department of Social Dentistry, USP/Araraquara
REGINA GUENKA PALMA-DIBB, Department of Operative Dentistry, USP/Ribeirão Preto



Coltene.com.br [/coltenebrasil](https://www.facebook.com/coltenebrasil) [Coltene Brasil](https://www.linkedin.com/company/coltene-brasil) [Coltene Brasil](https://www.youtube.com/coltenebrasil) [Coltene Brasil](https://www.instagram.com/coltenebrasil) [/coltenebrasil](https://twitter.com/coltenebrasil)



Vigodent S/A Indústria e Comércio
Rua Pesqueira, 26; CEP: 21041-150, Bonsucesso -
Rio de Janeiro / Brasil CNPJ 33.425.331/0001-22
sac@coltene.com.br - Hotline: 0800 0247422
MS Register No.: 10068870098
122012

COLTENE
Upgrade Dentistry