

WILLI GELLER  
*Creation*  
PRODUCT CATALOGUE





“I have no secret – my success is a result of all my experiences.”



Simply Brilliant – then as now

Inspired by the desire for a unique and vital copy of nature’s own creation, the internationally renowned dental technician Willi Geller created a metal and ceramics range in 1988 which made quite a splash: **perfect, almost natural restorations** became possible all at once – masterly in their physical properties and reliable and reproducible in terms of shading. With Creation, each tooth was suddenly able to retain its own character – in shade, form and function

The new ceramic materials from Switzerland quickly set the trend for dental technology worldwide. Geller’s **“master ceramics”** and his distinctive layering technique set benchmarks for a top-class dental aesthetic. Through the opportunity for ideal natural replicas, and through the seemingly **inexhaustible variety** of the metal and ceramics range, the dental technician gained an entirely new profile – the profile of an artist.

**Experience as a source of success**  
Founded as a “one-man company” 30 years ago, the dental company has continued to grow and develop steadily: new manufacturing technologies, additional fusions and effect materials, as well as **expanded indications and, more recently, a light-curing veneer composite with “Creation Brilliance”**. Today, Creation is one of the world’s most important manufacturers of ambitious ceramics systems, and resin and composite teeth with a strong profile. The innovative ceramic and composite materials and denture teeth are successfully processed in over 70 countries **across the globe**. And the users are not the only ones who are impressed; dentists and patients are also convinced by our “Simply Brilliant” aesthetics concept.

**Willi Geller’s dental creations represent the original – and are just as brilliant today.**

Creation Highlights

| VENEERING CERAMIC                                    | Page | CERAMIC ACCESSORIES    | Page |
|------------------------------------------------------|------|------------------------|------|
| <b>Metal</b>                                         |      | CreaPen-Separating Pen | 30   |
| Creation CC – Classic Metal Ceramic                  | 05   | Liquid Firing Pad      | 30   |
|                                                      |      | Creation Brush Kit     | 30   |
| <b>Zirconium</b>                                     |      | Aqualine master        | 31   |
| Creation ZI-CT – Zirconium Dioxide Ceramic           | 10   | Aqualine smart         | 31   |
| Creation ZI-F – Low-Fusing Zirconium Dioxide Ceramic | 14   | Aqualine mini          | 31   |
|                                                      |      |                        |      |
| <b>Lithium Disilicate</b>                            |      |                        |      |
| Creation LS/LS Press – Lithium Disilicate            |      |                        |      |
| Veneering and Pressable Ceramic                      | 19   |                        |      |
|                                                      |      |                        |      |
| COLOURS                                              | Page | COMPOSITE              | Page |
| <b>CreaColor</b>                                     |      | <b>Classic</b>         |      |
| Make up Neo                                          | 24   | Creation VC            | 33   |
| In Nova Neo                                          | 24   |                        |      |
|                                                      |      |                        |      |
| LIQUIDS                                              | Page | DENTURE TEETH          | Page |
| Modelling Liquid – High & Low Fusing                 | 26   | Creopal                | 41   |
| Opaque Liquid – High & Low Fusing                    | 26   | Creapearl 2            | 46   |
| Shoulder Liquid                                      | 26   |                        |      |
| Universal & Glaze Liquid                             | 27   |                        |      |
| Moisturing Fluid                                     | 27   |                        |      |
| In Nova Neo Fluid                                    | 27   |                        |      |
| Make up Neo Fluid                                    | 27   |                        |      |
| Make up Neo Reconditioner                            | 27   |                        |      |
|                                                      |      |                        |      |
| SHADE GUIDES                                         | Page |                        |      |
| Creation Shade Guides                                | 28   |                        |      |


# CERAMIC

“Creation is the only eligible choice for a dental technician to create convincing, which means naturally appearing, restorations.”

*Hubert Schenk*



REM-picture,  
etched, ceramic surface:



Leucite-reinforced metal ceramic CC

## CREATION CC – technically refined – inspired by nature

Creation CC is a high- fusing metal ceramic with unique optical and physical properties that has developed constantly in line with growing aesthetic demands – **recognized worldwide for 30 years!**

The perfectly coordinated ceramic materials have excellent homogeneity and thus guarantee high flexural strength. The result: a densely sintered structure for **pure and non-porous layering** with the utmost reliability.

The pure potash feldspar with micro-pure leucite crystal structure in Creation CC creates iridescent light refraction and natural brilliance and enhances the strength. Exceptional colour nuances can additionally be created with the unique effect materials.

The ceramic bonder Crea Alloy Bond also acts as a CTE buffer between metal and ceramic – for improved adhesive properties to non-precious metal and greater colour reliability.

With the modular classic brand Creation CC, metal-ceramic restorations can be produced with impressive aesthetics and dimensional stability – for **superb dental artistry!**

**Creation CC – natural brilliance and efficiency.**

## THE BRILLIANCE OF CREATION CC:

- Unique, lifelike aesthetics thanks to highly pure potash feldspar
- Resilience due to consistently high flexural strength
- Natural colour effect and light dynamics due to leucite crystals
- Ease of handling thanks to consistent colour and layering system
- Processing reliability since 1988



CC Creapast Kit | 000300  
*Complete paste opaques with modifier*

- 16 Opaque O-A1–D4 . . . . . 3g
- 5 Opaque Dentine Modifier . 3g
- 1 Creapast Brush
- 1 CC Manual



CC Bleach Shades Kit | 000250  
*Complete bleach shades colours*

- |                                               |                                   |
|-----------------------------------------------|-----------------------------------|
| 1 Bleach Creapast Opaque<br>O-AB . . . . . 3g | 1 Shoulder Liquid . . . . . 25 ml |
| 3 Bleach Dentine<br>BD-A, B, BO . . . . . 20g | 1 Modelling Liquid ML-H . 25 ml   |
| 1 Bleach Enamel S-AB . . . 20g                | 1 Bleach Shade Guide              |
| 1 Bleach Shoulder SP-AB . 15g                 | 1 CC Color Chart/Firing Chart     |



CC Opaque Kit | 000800  
*Complete powder opaques with modifier*

- |                                                |                  |
|------------------------------------------------|------------------|
| 16 Opaque O-A1–D4 . . . . . 20g                | 1 Shade Guide    |
| 8 Opaque Modifier<br>OM-11–17, G . . . . . 15g | 1 CC Color Chart |
| 1 Foundation Opaque<br>FOP . . . . . 20 g      | 1 CC Manual      |
| 2 Opaque Liquid OL-H . . 50ml                  |                  |
| 1 Universal & Glaze Liquid 25ml                |                  |



CC Gingiva Kit | 001590  
*Complete gingiva colours*

- |                                 |                                 |
|---------------------------------|---------------------------------|
| 1 Creapast Opaque O-G . . . 3g  | 1 Opaque Liquid OL-H . . 25 ml  |
| 1 Opaque Modifier OM-G . 15g    | 1 Modelling Liquid ML-H . 25 ml |
| 8 Gingiva G1–7, N . . . . . 15g | 1 Shade Guide                   |



CC Kit 1 | 000804  
*Basic porcelains kit (opaque dentines, dentines, incisals)*

- |                                             |                                  |
|---------------------------------------------|----------------------------------|
| 16 Opaque Dentine<br>OD-A1–D4 . . . . . 15g | 2 Modelling Liquid ML-H . 50 ml  |
| 16 Dentine D-A1–D4 . . . . . 20g            | 1 Universal & Glaze Liquid 25 ml |
| 4 Enamel E-57–60 . . . . . 20g              | 2 Shade Guides                   |
| 2 Transpa NT, OT . . . . . 20g              | 1 CC Color Chart                 |
| 1 Clear CL-O . . . . . 20g                  | 1 CC Manual                      |
| 1 Glaze GL . . . . . 10g                    |                                  |



CC Shoulder Kit | 001503  
*Complete shoulder porcelains*

- 10 Shoulder Powder  
SP-21–29, G . . . . . 15g
- 1 Shoulder Liquid . . . . . 25 ml
- 1 CC Color Chart/Firing Chart



CC Kit 2 | 001508  
*Effect porcelains kit*

- |                                                    |                                  |
|----------------------------------------------------|----------------------------------|
| 5 Opaque Dentine<br>OD-32, 37, 41, 43, 44. . . 15g | 4 Pearl Enamel PS-0–3. . . 20g   |
| 1 Dentine D-CW . . . . . 20g                       | 6 Neck Transpa HT-51–56 15g      |
| 2 Approximal Dentine<br>AD-1, 2. . . . . 15g       | 6 Make In MI-61–66 . . . . . 15g |
| 1 Clear UC . . . . . 20g                           | 1 Correction Powder KM. . . 20g  |
| 5 Transpa TI-1–5 . . . . . 20g                     | 2 Modelling Liquid ML-H 50ml     |
| 8 Effect Enamel<br>SI-01–06, SO-10, 11 . . . 20g   | 1 Universal & Glaze Liquid 25ml  |
|                                                    | 2 Shade Guides                   |
|                                                    | 1 CC Color Chart                 |
|                                                    | 1 CC Manual                      |



CC Transition Dentine Kit | 000801  
*Basic porcelains kit (transition dentines, incisals)*

- |                                                 |                                  |
|-------------------------------------------------|----------------------------------|
| 16 Opaque Dentine<br>OD-A1–D4 . . . . . 15g     | 2 Modelling Liquid ML-H . 50 ml  |
| 16 Transition Dentine<br>TD-A1–D4 . . . . . 20g | 1 Universal & Glaze Liquid 25 ml |
| 4 Enamel E-57–60 . . . . . 20g                  | 1 Shade Guide                    |
| 1 Clear CL-O . . . . . 20g                      | 1 CC Color Chart                 |
| 2 Transpa NT, OT . . . . . 20g                  | 1 CC Manual                      |
| 1 Glaze GL . . . . . 10g                        |                                  |



## CC Starter Kit | 000915

Introduction kit with tooth shades A2, A3, A3.5, B2, B3

|   |                 |      |   |                          |       |
|---|-----------------|------|---|--------------------------|-------|
| 2 | Crea Alloy Bond | 4 g  | 5 | Transpa NT, OT,          |       |
| 5 | Creapast Opaque | 3 g  |   | TI-1, 4, 5               | 20 g  |
| 5 | Opaque Dentine  | 15 g | 2 | Effect Enamel SI-02, 04  | 20 g  |
| 3 | Opaque Dentine  |      | 3 | Neck Transpa HT-52-54    | 15 g  |
|   | OD-32, 37, 43   | 15 g | 1 | Glaze GL                 | 10 g  |
| 5 | Dentine         | 20 g | 2 | Make In MI-63, 65        | 15 g  |
| 2 | Enamel E-58, 59 | 20 g | 2 | Modelling Liquid ML-H    | 25 ml |
| 1 | Clear CL-O      | 20 g | 2 | Universal & Glaze Liquid | 25 ml |
|   |                 |      | 1 | CC Manual                |       |



## CC Introduction Kit | 000920

Introduction kit with tooth shades A1, A2, A3, A3.5, B1, B2, C1, C2, C3, D2, D3

|                                  |                                   |
|----------------------------------|-----------------------------------|
| 11 Creapast Opaque. . . . . 3g   | 3 Effect Enamel                   |
| 3 Opaque Dentine                 | SI-02, 04, SO-10 . . . . . 20g    |
| OD-32, 37, 43 . . . . . 15g      | 1 Pearl Enamel PS-1 . . . . . 20g |
| 11 Dentine . . . . . 20g         | 3 Neck Transpa HT-52-54 15g       |
| 4 Enamel E-57-60 . . . . . 20g   | 1 Glaze GL . . . . . 10g          |
| 1 Clear CL-O . . . . . 20g       | 1 Modelling Liquid ML-H . 50ml    |
| 4 Transpa NT, TI-1, 2, 4 . . 20g | 1 Universal & Glaze Liquid 25ml   |
|                                  | 1 CC Color Chart & CC Manual      |

## Crea Alloy Bond | 000080

Crea Alloy Bond is a bonding agent which was specially developed for neutralising leaking metallic oxides. At the same time it acts as an efficient CTE buffer between porcelain and metal alloy.



- Safety owing to high flexural strength (70 MPa)
- No marginal discolorations – staining metallic oxides are neutralised
- Greater colour reliability because of the golden yellow undertone
- Compatible with all materials from Creation CC and LF


## WOP Opaque

The coating of opaque can be carried out at only 950 °C instead of 980 °C with the WOP paste. This procedure is useful with certain high-gold alloys and galvanised frameworks. After the WOP coating the second coating of opaque is applied.




- Firing temperature 950 °C

|                        |    |        |
|------------------------|----|--------|
| CC Creapast WOP sand   | 3g | 000322 |
| CC Creapast WOP sahara | 3g | 000323 |
| CC Creapast WOP gold   | 3g | 000324 |







REM-pictures  
etched, ceramic surfaces:



Glass ceramic



Leucite-reinforced metal ceramic CC



Leucite-reinforced  
zirconium dioxide ceramics ZI-CT

CREATION ZI-CT –  
the fine art of ceramics

Creation ZI-CT is an innovative feldspar ceramic containing leucite crystals which is matched to the CTE of all high-strength zirconium dioxide frameworks and structures – whether white, coloured or translucent.

The high proportion of pure potash feldspars creates light scatter that is comparable to that of natural teeth and a **brilliant colour effect**, while the micro-fine leucite crystals improve mechanical strength and thereby ensure a **durable ceramic structure**.

Fascinatingly fluorescent: the special bonding agent Frame Shade NT guarantees strong adhesion to the zirconium dioxide framework and also increases its fluorescence – for **powerful bond strength** without tensile stresses.

The modular Creation ZI-CT allows restorative work on zirconium dioxide with impressive colour and dimensional stability and with **brilliant light reflection and dynamics**.

High-fusing, high aesthetics, high stability: Creation ZI-CT.

*“Unobtrusively beautiful – a true Creation!”*  
Andreas Nolte, MDT

THE BRILLIANCE OF  
CREATION ZI-CT

- Natural colour effect and light dynamics due to leucite crystals
- Excellent bonding strength thanks to high firing temperature of 910 °C
- Extremely low shrinkage for fewer correction firings and more efficiency
- Reliability due to homogeneous crystalline structure and high resilience



ZI-CT Kit 1 | 620501  
Basic porcelains kit (dentines, incisals)

- |                                                    |                                 |
|----------------------------------------------------|---------------------------------|
| 1 Frame Shade NT . . . . . 4g                      | 1 Correction Powder KM . 20g    |
| 5 Opaque Dentine<br>OD-32, 37, 41, 43, 44. . . 15g | 1 Glaze GL . . . . . 10g        |
| 16 Dentine D-A1–D4. . . . . 20g                    | 1 Universal & Glaze Liquid 25ml |
| 4 Enamel E-57–60 . . . . . 20g                     | 2 Modelling Liquid ML-H . 50ml  |
| 1 Clear CL-O . . . . . 20g                         | 1 ZI-CT Color Chart & Manual    |
| 2 Transpa NT, OT . . . . . 20g                     |                                 |
| 2 Approximal Dentine<br>AD-1, 2. . . . . 15g       |                                 |



ZI-CT Kit 2 | 620502  
Effect porcelains kit

- |                                                  |                                 |
|--------------------------------------------------|---------------------------------|
| 1 Clear UC . . . . . 20g                         | 1 Universal & Glaze Liquid 25ml |
| 5 Transpa TI-1–5. . . . . 20g                    | 2 Modelling Liquid ML-H . 50ml  |
| 8 Effect Enamel<br>SI-01–06, SO-10, 11 . . . 20g | 1 ZI-CT Color Chart & Manual    |
| 4 Pearl Enamel PS-0–3. . . 20g                   |                                 |
| 6 Neck Transpa<br>HT-51–56. . . . . 15g          |                                 |
| 6 Make In MI-61–66 . . . . 15g                   |                                 |



ZI-CT Starter Kit | 620503  
Introduction kit with tooth shades A1, A2, A3, A3.5, B2, B3

- |                                                 |                                  |
|-------------------------------------------------|----------------------------------|
| 1 Frame Shade NT . . . . . 4g                   | 1 Pearl Enamel PS-0. . . . . 20g |
| 3 Opaque Dentine<br>OD-32, 37, 43 . . . . . 15g | 3 Neck Transpa HT-52–54 15g      |
| 6 Dentine . . . . . 20g                         | 3 Shoulder SP-21, 23, 28 . 15g   |
| 2 Enamel E-58, 59 . . . . . 20g                 | 2 Make In MI-63, 65 . . . . 15g  |
| 1 Clear CL-O . . . . . 20g                      | 1 Glaze GL . . . . . 10g         |
| 6 Transpa NT, OT<br>TI-1, 3–5 . . . . . 20g     | 2 Modelling Liquid ML-H . 25ml   |
| 2 Effect Enamel SI-02, 04 20g                   | 1 Shoulder Liquid . . . . . 25ml |
|                                                 | 1 Universal & Glaze Liquid 25ml  |
|                                                 | 1 ZI-CT Color Chart & Manual     |



ZI-CT Introduction Kit | 620506  
Introduction kit with tooth shades A1, A2, A3, A3.5, B1, B2, C1, C2, C3, D2, D3

- |                                                 |                                         |
|-------------------------------------------------|-----------------------------------------|
| 1 Frame Shade NT . . . . . 4g                   | 1 Pearl Enamel PS-0. . . . . 20g        |
| 3 Opaque Dentine<br>OD-32, 37, 43 . . . . . 15g | 3 Neck Transpa<br>HT-52–54. . . . . 15g |
| 11 Dentine . . . . . 20g                        | 1 Glaze GL . . . . . 10g                |
| 1 Clear CL-O . . . . . 20g                      | 2 Modelling Liquid ML-H . 50ml          |
| 4 Enamel E-57–60 . . . . . 20g                  | 1 Universal & Glaze Liquid 25ml         |
| 4 Transpa NT, TI-1, 3, 4 . . 20g                | 1 ZI-CT Color Chart & Manual            |
| 1 Effect Enamel SO-10 . . . 20g                 |                                         |



ZI-CT Bleach Shades Kit | 620505  
*Complete bleach shades colours*

- 3 Bleach Dentine  
BD-A, B, BO . . . . . 20g
- 1 Bleach Enamel S-AB . . . 20g
- 1 Shoulder SP-AB. . . . . 15g
- 1 Modelling Liquid ML-H 25ml
- 1 Shoulder Liquid . . . . . 25ml
- 1 ZI-CT Color Chart




ZI-CT Gingiva Kit | 620504  
*Complete gingiva colours*

- 8 Gingiva G1–7, N. . . . . 15g
- 1 Modelling Liquid ML-H 25ml
- 1 ZI-CT Color Chart  
Bleach & Gingiva



ZI-CT Shoulder Kit | 620507  
*Complete shoulder porcelains*

- 10 Shoulder  
SP-21–29, G . . . . . 15g
- 1 Shoulder Liquid . . . . . 25ml
- 1 ZI-CT Color Chart/Firing Chart





CREATION ZI-F – naturally brilliant – brilliantly natural

Creation ZI-F is a tried and tested zirconium dioxide ceramic. It is perfectly suited to veneering any kind of zirconium dioxide framework. At the same time, the modular materials are also an ideal choice as **add-on materials** for Creation ZI-CT because their firing temperature is about 100 °C lower than the high-fusing zirconium dioxide veneering ceramic.

Creation ZI-F is an **extremely versatile** zirconium dioxide ceramic with high **natural light reflection**, combined with all the material-related advantages of zirconium dioxide. With the unique effect materials, exceptional colour nuances can also be created rapidly and comfortably – for maximum impact with minimum effort.

Creation ZI-F – natural light dynamics!



ZI-F Kit 1 | 653010  
Basic porcelains kit (dentines, incisals)

- |                                           |                                 |
|-------------------------------------------|---------------------------------|
| 5 Frame Shade FS-1-4, IL . . . . . 10g    | 2 Modelling Liquid ML-L . 50ml  |
| 16 Dentine D-A1–D4 . . . . . 20g          | 1 Universal & Glaze Liquid 25ml |
| 4 Enamel E-57–60 . . . . . 20g            | 2 Shade Guides                  |
| 2 Clear CL-O, UC . . . . . 20g            | 1 ZI-F Color Chart              |
| 6 Transpa NT, OT, TI-1, 3–5 . . . . . 20g | 1 ZI-F Manual                   |
| 1 Glaze GL . . . . . 10g                  |                                 |



ZI-F Kit 2 | 653015  
Effect porcelains kit

- |                                                |                                            |
|------------------------------------------------|--------------------------------------------|
| 5 Opaque Dentine OD-32, 37, 41, 43, 44 . . 15g | 2 Approximal Dentine AD-1, 2 . . . . . 15g |
| 5 Effect Enamel SI-02, 04, 06, SO-10, 11 20g   | 1 Shoulder Liquid . . . . . 25ml           |
| 2 Pearl Enamel PS-0, 3 . . . 20g               | 2 Modelling Liquid ML-L . 25ml             |
| 3 Neck Transpa HT-52–54 15g                    | 1 Universal & Glaze Liquid 25ml            |
| 5 Shoulder SP-22–25, 28 . 15g                  | 2 Shade Guides                             |
| 3 Make In MI-61, 63, 65 . . 15g                | 1 ZI-F Color Chart                         |
|                                                | 1 ZI-F Manual                              |
|                                                | 1 ZI-F Color Chart/Firing Chart            |



ZI-F Bleach Shades Kit | 652050  
Complete bleach shades colours

- |                                            |                                  |
|--------------------------------------------|----------------------------------|
| 1 Bleach Frame Shade FS-AB . . . . . 10g   | 1 Shoulder Liquid . . . . . 25ml |
| 3 Bleach Dentine BD-A, B, BO . . . . . 20g | 1 Modelling Liquid ML-L . 25ml   |
| 1 Bleach Enamel S-AB . . . 20g             | 1 Shade Guide                    |
| 1 Bleach Shoulder SP-AB . 15g              | 1 ZI-F Color Chart/Firing Chart  |



ZI-F Gingiva Kit | 603030  
Complete gingiva colours

- |                                 |
|---------------------------------|
| 1 Frame Shade FS-G . . . . 10g  |
| 8 Gingiva G1–7, N . . . . . 15g |
| 1 Modelling Liquid ML-L . 25ml  |
| 1 Shade Guide                   |



ZI-F Starter Kit | 653025  
Introduction kit with tooth shades A1, A2, A3, A3.5, B2, B3

- |                                              |                                  |
|----------------------------------------------|----------------------------------|
| 2 Frame Shade FS-1, 2 . . . 10g              | 3 Neck Transpa HT-52–54 15g      |
| 3 Opaque Dentine OD-32, 37, 43 . . . . . 15g | 3 Shoulder SP-22, 23, 28 . 15g   |
| 6 Dentine . . . . . 20g                      | 2 Make In MI-63, 65 . . . . 15g  |
| 2 Enamel E-58, 59 . . . . . 20g              | 1 Glaze GL . . . . . 10g         |
| 1 Clear CL-O . . . . . 20g                   | 1 Shoulder Liquid . . . . . 25ml |
| 6 Transpa NT, OT, TI-1, 3–5 20g              | 2 Modelling Liquid ML-L . 25ml   |
| 2 Effect Enamel SI-02, 04 . 20g              | 1 Universal & Glaze Liquid 25ml  |
| 1 Pearl Enamel PS-02 . . . 20g               | 1 ZI-F Color Chart               |
|                                              | 1 ZI-F Manual                    |



ZI-F Single Products

- |                       |                      |
|-----------------------|----------------------|
| Frameshade FS-NT      | 601087 . . . . . 4g  |
| Shoulder Powder SP-21 | 602121 . . . . . 15g |
| Transpa TI-2          | 652077 . . . . . 20g |
| Neck Transpa HT-51    | 652151 . . . . . 15g |
| Neck Transpa HT-56    | 652156 . . . . . 15g |
| Correction Powder KM  | 652090 . . . . . 20g |

16 Creation ZI-F – Color Chart

A1

A2

A3

A3,5

A4

B1

B2

B3

B4

C1

C2

C3

C4

D2

D3

D4

DENTINE

OD-32  
havanna

OD-37  
curry

OD-41  
orange

OD-43  
ivory

OD-44  
cuba

57

58

59

60

OPAQUE DENTINEENAMEL

CL-O  
clear

UC  
window

NT  
neutral

OT  
opal

TI-01  
blue

TI-02  
white

TI-03  
pink

TI-04  
yellow

TI-05  
grey

CLEARTRANSPA

PS-0  
white

PS-3  
orange

SI-02  
medium  
yellow

SI-04  
light  
orange

SI-06  
heavy  
orange

SO-10  
blue

SO-11  
orange

HT-51  
irides-  
cent

HT-52  
khaki

HT-53  
sand

HT-54  
honey

HT-56  
ocher

PEARL ENAMELEFFECT ENAMELNECK TRANSPA HIGH FLUORESCENT

SP-21  
neutral

SP-22  
flamingo

SP-23  
sand

SP-24  
gold

SP-25  
red  
brown

SP-28  
opaque  
yellow

SHOULDER POWDERSHIGH FLUORESCENT

MI-61  
ivory

MI-63  
honey  
yellow

MI-65  
gold

AD-1  
light  
yellow

AD-2  
orange

G1  
purple

G2  
dark  
pink

G3  
light  
pink

G4  
flamingo

G5  
rose

G6  
dark pink  
opaque

G7  
violet

G-N  
neutral

MAKE IN  
HIGH FLUORESCENT  
INTERNAL POWDERAPPROXIMAL D.GINGIVAL

BD-A

BD-B

BD-BO

S-AB

SP-AB

KM

BLEACH DENTINEBLEACH EN.BLEACH SH.CORRECTION POWDER

Shades A1–D4 are referring to the Vita® Shade Guide. Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.



# LITHIUM DISILICATE

“With Creation LS, there is now finally a veneering ceramic for lithium disilicate, which offers me a wide scope for creating restorations with a natural effect. The colour concept matches that of the tried and tested Creation ceramics, with an outstanding light reflection, plenty of brightness and strong, dynamic colours. The combination of Creation LS Press and Creation LS provides me with many benefits in my daily work, and makes it easier for me to create restorations with a vivid effect.”

*Sebastien Mosconi, Oral Design*



## CREATION LS/LS Press – Press, veneer and individualise with lithium disilicate!

Efficiency and natural aesthetics: Creation LS Press and Creation LS are a perfect ceramic symbiosis offering dental technicians reliability and stability but also brilliant colours and versatility on and with lithium disilicate, as standard is not enough for Creation!

Creation Willi Geller completes the product range with a new pressable and veneering ceramic allowing the dental technician outstanding aesthetics, whether used single or in combination. With the Creation LS Press blanks and the Creation LS veneering ceramic perfectly matched to the lithium disilicate frame highly aesthetic and stable all-ceramic restorations can be manufactured – efficiently pressed and then individually veneered resp. characterised with the proven CreaColor stains and modifiers.

Reliably combined for extended possibilities: the two system components Creation LS Press and Creation LS follow master ceramist Willi Geller's "Simply Brilliant" philosophy proven for 30 years. They allow the dental technician manufacturing naturally beautiful anterior and posterior restorations with resp. on lithium disilicate assuring system safety.

## Creation LS / LS Press – the sophisticated combination

### THE BRILLIANCE OF THE CREATION LS:

- Comprehensive colour assortment for individual veneering and characterisation
- Layer build-up according to the tried-and-tested colour coding concept of Creation Willi Geller
- Natural colour brilliance and light effects due to feldspar content
- Homogeneous crystal structure and outstanding surface density
- Maximum stability – even after several firings

### THE BRILLIANCE OF THE CREATION LS PRESS:

- High strength due to homogeneous crystalline structure: > 500 MPa
- Minimal invasive preparation: 0.3 mm veneer, 0.8 mm anterior
- Creative freedom of manufacturing: full veneers, partial veneers
- The appropriate blank for every indication



- LS Kit | 830083**  
*Basic and effect porcelains with the tooth shades BD-A, BD-B, BD-B0, A1, A2, A3, A3.5, B1, B2, B3, C2, D2*
- |                                                              |                                                    |
|--------------------------------------------------------------|----------------------------------------------------|
| 4 Opaque Dentine<br>OD-32, 37, 41, 43 . . . . . 15 g         | 3 Effect Enamel<br>SI-02, 04, SO-10 . . . . . 20 g |
| 12 Dentine BD-A, BD-B, BD-B0,<br>A1-A3.5, B1-B3, C2, D2 20 g | 1 Pearl Enamel PS-0 . . . . . 20 g                 |
| 4 Enamel E-57-60 . . . . . 20 g                              | 4 Neck Transpa HT-51-54. 15 g                      |
| 1 Bleach Enamel S-AB . . . 20 g                              | 2 Make-in MI-61, 63 . . . . . 15 g                 |
| 1 Clear CL-0 . . . . . 20 g                                  | 1 Glaze GL . . . . . 10 g                          |
| 7 Transpa NT, OT,<br>TI-01-05 . . . . . 20 g                 | 2 Modelling Liquid ML-L . 50 ml                    |
|                                                              | 1 Universal & Glaze Liquid 25 ml                   |
|                                                              | 1 Creation LS/LS Press Manual                      |



- LS Press Starter Kit | 830080**  
*Introduction kit with the tooth shades BD-A, BD-B0, A1, A2, A3 and matching press blanks with different opacities MO, LT, MT, HT*
- |                                                  |                                                |
|--------------------------------------------------|------------------------------------------------|
| 3 LS Press MO-0-2. . . . 5 x 3 g                 | 1 Clear CL-0 . . . . . 20 g                    |
| 1 LS Press LT-A. . . . . 5 x 3 g                 | 1 Transpa OT . . . . . 20 g                    |
| 1 LS Press MT-B00. . . . 5 x 3 g                 | 2 Effect Enamel<br>SI-02, SO-10 . . . . . 20 g |
| 2 LS Press HT-2, 3 . . . . 5 x 3 g               | 1 Pearl Enamel PS-0 . . . . 20 g               |
| 2 Opaque Dentine<br>OD-32, 41 . . . . . 15 g     | 1 Make-in MI-61, 63 . . . . 15 g               |
| 5 Dentine BD-A, BD-B0,<br>D-A1-A3 . . . . . 20 g | 1 Glaze GL . . . . . 10 g                      |
| 2 Enamel E-58, 59 . . . . . 20 g                 | 1 Modelling Liquid ML-L . 50 ml                |
| 1 Bleach Enamel S-AB . . . 20 g                  | 1 Universal & Glaze Liquid 25 ml               |
|                                                  | 1 Creation LS/LS Press Manual                  |



- LS Press Introduction Kit A1, A2 | 830081**  
*The easy introduction into the system with the tooth shades A1, A2 and the matching press blanks with different opacities MO-1, LT-A, HT-2*
- |                                   |                                  |
|-----------------------------------|----------------------------------|
| 1 LS Press MO-1 . . . . . 5 x 3 g | 1 Clear CL-0 . . . . . 20 g      |
| 1 LS Press LT-A. . . . . 5 x 3 g  | 1 Transpa OT . . . . . 20 g      |
| 1 LS Press HT-2 . . . . . 5 x 3 g | 1 Glaze GL . . . . . 10 g        |
| 1 Opaque Dentine OD-32 . 15 g     | 1 Modelling Liquid ML-L . 50 ml  |
| 2 Dentine D-A1, A2 . . . . . 20 g | 1 Universal & Glaze Liquid 25 ml |
| 1 Enamel E-58 . . . . . 20 g      | 1 Creation LS/LS Press Manual    |



- LS Press Introduction Kit A3 | 830082**  
*The easy introduction into the system with the tooth shade A3 and the matching press blanks with different opacities MO-2, LT-A, HT-3*
- |                                   |                                  |
|-----------------------------------|----------------------------------|
| 1 LS Press MO-2 . . . . . 5 x 3 g | 1 Clear CL-0 . . . . . 20 g      |
| 1 LS Press LT-A. . . . . 5 x 3 g  | 1 Transpa OT . . . . . 20 g      |
| 1 LS Press HT-3 . . . . . 5 x 3 g | 1 Glaze GL . . . . . 10 g        |
| 1 Opaque Dentine OD-32 . 15 g     | 1 Modelling Liquid ML-L . 50 ml  |
| 2 Dentine D-A3 . . . . . 20 g     | 1 Universal & Glaze Liquid 25 ml |
| 1 Enamel E-59 . . . . . 20 g      | 1 Creation LS/LS Press Manual    |
| 1 Effect Enamel SO-10 . . . 20 g  |                                  |



**Creation LS Press MO-Ingot (Medium Opacity)**  
MO press ingots are the first choice for full veneering. Due to their high opacity, they feature a high light reflection and are suitable for covering discoloured dies and titanium abutments.

**Colours**  
MO-0, MO-1, MO-2 . . . . . 5 x 3 g



**Creation LS Press LT-Ingot (Low Translucency)**  
LT press ingots can be used for restorations to be fully or partially veneered with Creation LS. With thin layer thicknesses, Creation LT can also be used for monolithic restorations.

**Colours**  
LT-A, LT-B, LT-C, LT-D . . . . . 5 x 3 g



**Creation LS Press MT-Ingot (Medium Translucency)**  
MT press ingots are the first choice for fully-anatomic restorations which are characterised by individual staining. The MT ingot can also be used for full or partial veneering.

**Colours**  
MT-B00, MT-B0, MT-A1-3, MT-B1, B2, MT-D2 . . . . . 5 x 3 g



**Creation LS Press HT-Ingot (High Translucency)**  
HT press ingots are particularly suitable for minimal invasive restorations (e.g. inlays, onlays and veneers) where mainly the natural enamel is replaced. The shade individualisation is carried out in the staining technique.


**Colours**  
HT-B00, HT-B0, HT-1-4 . . . . . 5 x 3 g

Processing techniques


| Opacity                   | Indications      |         |        |        |       | Processing techniques |                    |                    |
|---------------------------|------------------|---------|--------|--------|-------|-----------------------|--------------------|--------------------|
|                           | Occlusal veneers | Veneers | Inlays | Onlays | Crown | Staining technique    | Cut-back technique | Layering technique |
| MO<br>Medium Opacity      |                  |         |        |        | ✓     |                       |                    | ✓                  |
| LT<br>Low Translucency    |                  |         |        |        | ✓     |                       | ✓                  | ✓                  |
| MT<br>Medium Translucency | ✓                | ✓       | ✓      | ✓      | ✓     | ✓                     | ✓                  |                    |
| HT<br>High Translucency   | ✓                | ✓       | ✓      | ✓      |       | ✓                     |                    |                    |




Layering technique



Cut-back technique



Staining technique




Colour reference table

| Tooth colour              | A1    | A2    | A3        | A3,5 | A4    | B1        | B2    | B3    | B4    | C1    | C2    | C3    | C4    | D2    | D3    | D4    |
|---------------------------|-------|-------|-----------|------|-------|-----------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| MO<br>Medium Opacity      | MO-1  |       | MO-2      |      |       | MO-1      |       | MO-2  |       | MO-1  |       | MO-2  |       |       |       |       |
| LT<br>Low Translucency    | LT-A  |       |           |      |       | LT-B      |       |       |       | LT-C  |       |       | LT-D  |       |       |       |
| MT<br>Medium Translucency | MT-A1 | MT-A2 | MT-A2/-A3 |      | MT-A3 | MT-B0/-B1 | MT-B2 | MT-B2 | MT-A3 | MT-B1 | MT-A1 | MT-D2 | MT-D2 | MT-D2 | MT-D2 | MT-A2 |
| HT<br>High Translucency   | HT-2  |       | HT-3      |      | HT-4  | HT-1      | HT-3  |       |       | HT-4  | HT-3  |       | HT-4  |       | HT-3  |       |

Materials which are not related to the Vita®shade guide:

|        |      |        |       |        |       |
|--------|------|--------|-------|--------|-------|
| Bleach | MO-0 | MT-B00 | MT-B0 | HT-B00 | HT-B0 |
|--------|------|--------|-------|--------|-------|





**CREACOLOR –**  
stains and modifiers for all ceramics

CreaColor is a complete **“two in one” colour system** which can be used to place individual colour nuances simply and efficiently **for all ceramic surfaces and materials** – whether by the internal or external staining technique. The physical properties of the **high-quality stains and modifiers** are so perfectly matched to the ceramic surfaces and materials that they can be worked universally with all layering and pressed ceramics, pressed or fired dentine cores as well as zirconium dioxide frames – after all, **colour creativity should be unlimited!** Furthermore, the new stains are extremely fine-grained and fluorescent and they can readily be mixed with each other.

CreaColor is made up of two components:

- CreaColor Make up Neo
- CreaColor In Nova Neo

**THE BRILLIANCE OF  
CREACOLOR:**

- **Utmost colour brilliance for internal and external colouring**
- **Versatile and efficient: one system for all ceramics**
- **Natural colour effect due to very high fluorescence values**
- **Easy to handle: can be mixed with each other and can be corrected**



CreaColor Kit | 001200  
Contains Make up Neo Kit & In Nova Neo Kit



Make up Neo Kit | 001100  
Glazing and fluorescent stains

- 13 Fluorescent Stain Glazes  
MUN-F1–F13 . . . . . 4g
- 4 Dentine Shades MUN-FA–FD  
in the four Vita® colour groups for colour corrections  
and for painting onto press ingots . . . . . 4g
- 1 Fluorescent Glaze GL-F . . . . . 4g
- 1 Make up Neo Fluid . . . . . 25ml
- 1 Make up Neo Reconditioner Fluid . . . . . 8ml  
(for freshening up concentrated stains)
- 1 Brush No. IV



In Nova Neo Kit | 001150  
All-purpose modifiers

- 10 Fluorescent Modifiers  
INN-1–10. . . . . 3g
- 2 Effect Shades IL, CRL (Illusion and Crackliner)  
for individual accents . . . . . 3g
- 1 In Nova Neo Fluid . . . . . 25ml
- 1 Brush No. IV

CREACOLOR COLOR CHART

MAKE UP NEO

MUN-F1  
white

MUN-F2  
eggshell

MUN-F3  
lemon  
yellow

MUN-F4  
mandarine

MUN-F5  
flamingo

MUN-F6  
apricot

MUN-F7  
light brown

MUN-F8  
medium  
brown

MUN-F9  
olive

MUN-F10  
red

MUN-F11  
blue

MUN-F12  
grey

MUN-F13  
deep blue

MUN-FA  
stain A

MUN-FB  
stain B

MUN-FC  
stain C

MUN-FD  
stain D

MUN-GL-F  
Glaze

UF

Recon-  
ditionner

IN NOVA NEO

INN-1  
melone

INN-2  
rose

INN-3  
yellow

INN-4  
red  
brown

INN-5  
brown

INN-6  
olive

INN-7  
grey

INN-8  
blue

INN-9  
white

INN-10  
mais

INN-CRL  
Crackli-  
ner

INN-IL  
Illusion

INN-  
FLUID

**Fewer bottles – more simplicity and efficiency:** For more clarity and space on the laboratory table, Creation liquids can be obtained in high fusion or low fusion versions.



**Modelling Liquid – High & Low Fusing**  
Standard modelling liquid for mixing Creation ceramics except for shoulder ceramics (SP).

- ML-H (High) for high fusion ceramics Creation CC and ZI-CT
- ML-L (Low) for low fusion ceramic Creation ZI-F

Modelling Liquid High Fusing . . . . . 25 ml | 000012  
Modelling Liquid High Fusing . . . . . 50 ml | 000013  
Modelling Liquid High Fusing . . . . . 250 ml | 000014  
Modelling Liquid High Fusing . . . . . 1.000 ml | 000015

Modelling Liquid Low Fusing . . . . . 25 ml | 000016  
Modelling Liquid Low Fusing . . . . . 50 ml | 000017  
Modelling Liquid Low Fusing . . . . . 250 ml | 000018



**Opaque Liquid – High Fusing**  
Mixing liquid for the high fusion ceramic Creation CC.

Opaque Liquid High Fusing . . . . . 25 ml | 000019  
Opaque Liquid High Fusing . . . . . 50 ml | 000020  
Opaque Liquid High Fusing . . . . . 250 ml | 000021



**Shoulder Liquid**  
Universal modelling liquid for mixing all SP shoulder ceramics from Creation.

Shoulder Liquid . . . . . 25 ml | 000025  
Shoulder Liquid . . . . . 50 ml | 000026



**Universal & Glaze Liquid**  
Additional universal liquid with a wide range of uses.

- Addition of a very small amount into the ceramic mixed with ML enhances the smoothness
- The pre-drying time must be extended during firing
- Mixing all Creation glazing materials and CreaColor colours In Nova Neo and Make up Neo
- Remoistening of Creapast

Universal & Glaze Liquid . . . . . 25 ml | 000010  
Universal & Glaze Liquid . . . . . 50 ml | 000011



**Moisturing Fluid**  
For extensive restorations and complex corrections, as a moisturiser for the applied ceramics.

- No more drying out, the ceramics remain consistently modellable
- Efficient and simple: one fluid for all ceramics
- Longer working time for mixed materials – especially for sizeable work
- Targeted application due to controllable moisture content

Moisturing Fluid . . . . . 25 ml | 000090



**In Nova Neo Fluid**  
Mixing liquid for In Nova Neo.

In Nova Neo Fluid . . . . . 25 ml | 001163



**Make up Neo Fluid**  
Mixing liquid for Make up Neo.

Make up Neo Fluid . . . . . 25 ml | 001125

**Make up Neo Reconditioner**  
For freshening up concentrated stains.

Make up Neo Reconditioner . . . . . 8 ml | 001126

## KIT-BASED SHADE GUIDES

| Art. no. | Article name                             |
|----------|------------------------------------------|
| 000064   | CC/Shade Guide Kit 1 804                 |
| 000073   | CC/Shade Guide Kit 2 1508                |
| 000076   | ZI-F.LF.TI/Shade Guide Kit 1 (28 Shades) |
| 000077   | ZI-F.LF.TI/Shade Guide Kit 2 (26 Shades) |
| 000087   | CC/Shade Guide Transition Dentine        |
| 000092   | ZI-CT/Shade Guide Kit 1                  |
| 000093   | ZI-CT/Shade Guide Kit 2                  |

## SHADE GUIDES

Shades A1–D4 are referring to the Vita® Shade Guide.  
Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

**Shade Guide**  
**Dentine A-D**

Art. no. 000040

Shade Guide  
**Opaque Dentine**  
**A-D**


Art. no. 000041



|    |    |    |      |    |    |    |    |    |
|----|----|----|------|----|----|----|----|----|
| A1 | A2 | A3 | A3,5 | A4 | B1 | B2 | B3 | B4 |
| C1 | C2 | C3 | C4   | D2 | D3 | D4 |    |    |

Shade Guide  
**Transpa**

Art. no. 000042



The diagram shows a cross-section of a dental crown with five shade tabs on its top surface. The tabs are labeled CL-O, UC, NT, OT, and KM from left to right. The crown has a central circular opening and a surrounding ring.

Shade Guide  
**Opaque Dentine  
+ AD**

Art. no. 000043



The shade guide consists of seven circular color swatches arranged horizontally. From left to right, the colors are: a light yellow, a slightly darker yellow, a light orange, a very light cream, a light brown, a light yellow, and a light yellow. Below each swatch is a label: OD-32, OD-37, OD-41, OD-43, OD-44, AD-1, and AD-2.

| OD-32 | OD-37 | OD-41 | OD-43 | OD-44 | AD-1 | AD-2 |
|-------|-------|-------|-------|-------|------|------|
|       |       |       |       |       |      |      |

**Shade Guide**  
**Enamel**

Art. no. 000044


|       |       |       |      |       |       |       |       |       |
|-------|-------|-------|------|-------|-------|-------|-------|-------|
| 57    | 58    | 59    | 60   | SI-01 | SI-02 | SI-03 | SI-04 | SI-05 |
| SI-06 | SO-10 | SO-11 | PS-0 | PS-1  | PS-2  | PS-3  |       |       |

Shade Guide  
**Transpa Modifier**

Art. no. 000045

|       |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| TI-01 | TI-02 | TI-03 | TI-04 | TI-05 |       |
| HT-51 | HT-52 | HT-53 | HT-54 | HT-55 | HT-56 |

Shade Guide  
**Make In**



Art. no. 000046

**Shade Guide**

**Shoulder**

Art. no. 000047

|       |       |       |       |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| SP-21 | SP-22 | SP-23 | SP-24 | SP-25 | SP-26 | SP-27 | SP-28 | SP-29 |
|       |       |       |       |       |       |       |       |       |
|       |       |       |       |       |       |       |       |       |

Shade Guide  
**Gingiva**



Art. no. 000048



| G1 | G2 | G3 | G4 | G5 | G6 | G7 | G-N |
|----|----|----|----|----|----|----|-----|
|    |    |    |    |    |    |    |     |

Shade Guide  
**Bleach**

Art. no. 000049

BD-A BD-B BD-BO S-AB SP-AB

Shade Guide  
**Wheel Empty**

Art. no. 000091

Shade Guide  
**Single Shade Tab**



Art. no. 000058 + Color description



**CreaPen Separating Pen**  
The double-ended ceramic insulation pen with one pointed and one broad nib for surfaces or grooves respectively can be simply used on the sealed die.

CreaPen Separating Pen | 000059



**Liquid Firing Pad**  
The liquid firing pad for all Creation systems precisely fixes the object on the fireproof tray.

- Stable hold on the firing tray
- Easy to handle because of proven dispenser design
- Safe because of firing stability and optimum heat/warmth build-up onto the object
- Economical because of precise measuring
- Can be used for all Creation systems
- Tough and easy to remove

Liquid Firing Pad | 000081



**Creation Brush Kit**  
The six different sizes of the first-class Kolinsky red sable brush with a heavy wooden holder perfectly meet all the requirements of layering techniques.

- Creation-Brush-Kit | 000136
- |                          |                                  |
|--------------------------|----------------------------------|
| 1 Brush No. I   000131   | 4 Brush No. IV   000134          |
| 2 Brush No. II   000132  | 5 Brush No. V   000135           |
| 3 Brush No. III   000133 | 6 Brush Creapast No. VI   000321 |



**aqua line**

**Aqualine Mixing Trays**  
With Aqualine, materials will always have the right viscosity. Drying out is prevented by the cellulose strips continuously absorbing liquid from the reservoir.

**Aqualine master**  
Dimension: 335 x 225 x 25 mm  
Components: Tray with glass plate (light grey) and hinged cover, self-adhesive write-on labels, cellulose strips

Aqualine master . . . . . anthracite | AMA-001  
Aqualine master . . . . . white | AMA-002



**Aqualine smart**  
Dimension: 260 x 170 x 32 mm  
Components: Tray with glass plate (light grey) and hinged cover, self-adhesive write-on labels, cellulose strips

Aqualine smart . . . . . anthracite | AS-001  
Aqualine smart . . . . . white | AS-002

**Aqualine mini**  
Dimension: 133 x 71 x 23 mm  
Components: Tray with glass plate (light grey) and hinged cover, cellulose strips

Aqualine mini . . . . . anthracite | AM-001

ACCESSORIES

**Aqualine master** Glass Plate, printed light grey. . . . . 1 piece | AP-010

**Aqualine master** Cellulose Strips 40 pieces . . . . . AP-013

**Aqualine master** Labels, self-adhesive for stain . . . . . 5 pieces | AMA-018 for ceramic . . . . . 30 pieces | AMA-019

**Aqualine master/smart** Stain Kit, made of porcelain 1 stain tray/1 mixing tray. . . . . AS-016

**Aqualine smart** Glass Plate, printed light grey. . . . . 1 piece | AS-010

**Aqualine smart** Cellulose Strips 30 pieces . . . . . AS-012

**Aqualine mini** Tray Box with stain tray, with 15 cavities 1 piece . . . . . AM-002

**Aqualine mini** Cellulose Strips 40 pieces . . . . . AM-010

**Aqualine** Pool Set Dimension: 105 x 65 x 48 mm, rinsing cup 250 ml, sponge, sponge holder transparent. . . . . AS-014

**Aqualine** Pool Set Sponge 5 pieces . . . . . AS-015

**Aqualine** Glass Care 2 x 30 ml . . . . . AC-001

**Aqualine** Two in One Build up Liquid 500 ml . . . . . AC-002

# COMPOSITE

“Creation VC is an aesthetic combination offering us versatile solutions. At the same time it gives us complete freedom to create amazing results with relatively little effort, for all indications from the aesthetic temporary to the final work.”

*Hans-Joachim Lotz*



## CREATION VC – definitely aesthetic and versatile!

**Brilliant colours, natural light dynamic**, and high processing and system quality. The new micro-hybrid composite not only possesses unmatched durability and mechanical stability, but also well-balanced opalescence and translucence. Based on **hybrid MFR technology** and the **tried-and-tested colour coding concept**, it offers dental technicians a unique, versatile, colourfast and sturdy composite system, with a total of 102 individual porcelains in a user-friendly packaging concept. The **flow materials** are a particular highlight and can be processed almost like ceramic materials due to their excellent application and modelling properties, as well as the wide range of natural gingival shades and make-up stains – for vivid and effective **aesthetics in red and white!**

On the basis of a unique filler structure and the familiar colour and layer system used by ceramics specialists, natural analogue composite restoration with the tried-and-tested Creation brilliance is now possible.

Thanks to the non-adhesive consistency and balanced viscosity, the new composite material is ideal for modelling and can be polished to a fine shine.

The modular **Creation VC** complete system offers the dental technician a virtually unlimited number of colour and structural combinations, a wide spectrum of different consistencies and processing possibilities, and a customised portfolio of new technologies such as cuvette technology (pressing and injection methods), as well as compatibility with **Creopal Shell**.

**Creation VC – an expressive and high-performance composite system.**

## THE BRILLIANCE OF CREATION VC:

- **Unique, natural-looking aesthetics** due to microhybrid filler structure/hybrid MFR technology
- **A wide variety of indications** for metal-supported and metal-free indications
- **Perfect, safe and co-ordinated combinations** with **Creopal Shells**
- **High mechanical strength**
- **Resilience** through equally high resistance to bending
- **Different viscosities**
- **Natural colour effect and light dynamics**
- **Safe processing and bonding** thanks to latest generation primers

VC Kit 1 | 820201  
Basic materials kit (opaques, opaque dentines, dentines, incisals)

**Amazingly simple and varied:**  
Beside 16 opaques, 5 opaque dentines and 16 standard dentines in high density colours, dental technicians also have 4 assigned incisal materials at their disposal for an individual aesthetic. In addition, the set contains a further 2 primers, a variety of accessories and 2 opaques:  
**VC Base Opaque O-BO:** Paste opaque with excellent flow and light-curing qualities, also sets in undercuts for mechanical retentions, and provide a good primer.  
**VC Cervical Opaque O-CO:** effectively prevents metal edges from showing through, directly applied after the VC Base Opaque. VC Cervical can also be used occlusally or lingually for intensifying, or mixed with the opaques.



- |                                   |                              |
|-----------------------------------|------------------------------|
| 1 Base Opaque O-BO. . . 2.4 ml    | 1 Metal Primer. . . . . 5 ml |
| 1 Cervical Opaque O-CO. 2.4 ml    | 10 Disposable Palette        |
| 16 Opaque O-A1–D4 . . . . 2.4 ml  | 1 Mixing Pad                 |
| 5 Opaque Dentine                  | 10 Brush Round, Size 1       |
| OD-32, 37, 41, 43, 44. . 2.9 ml   | 10 Brush Flat, Size 1        |
| 16 Dentine D-A1–D4 . . . . 2.9 ml | 2 Brush Holder               |
| 4 Enamel E-57–60 . . . . . 2.9 ml | 1 Brush No. 7                |
| 1 Separator. . . . . 5 ml         | 2 Syringe Plunger            |
| 1 Diamond Polisher. . . . . 8 g   | 3 Instructions for use       |
| 1 Air Protect . . . . . 10 ml     | 1 Creation VC Workbook       |
| 1 Composite Primer . . . . . 3 ml |                              |

VC Kit 1 Starter | 820206  
Introduction kit in the 5 most common base shades A2, A3, A3.5, B2, B3

Besides 5 opaques, 3 opaque dentines and 5 standard dentines in high density colours, dental technicians also have 3 assigned incisal materials at their disposal for individual aesthetics. In addition, the set contains a variety of accessories and the Composite Primer.



- |                                   |                                   |
|-----------------------------------|-----------------------------------|
| 1 Base Opaque O-BO. . . 2.4 ml    | 1 Air Protect . . . . . 10 ml     |
| 1 Cervical Opaque O-CO. 2.4 ml    | 1 Composite Primer . . . . . 3 ml |
| 5 Opaque                          | 10 Disposable Palette             |
| O-A2–A3.5, B2, B3 . . . 2.4 ml    | 1 Mixing Pad                      |
| 3 Opaque Dentine                  | 10 Brush Round, Size 1            |
| OD-32, 37, 43 . . . . . 2.9 ml    | 10 Brush Flat, Size 1             |
| 5 Dentine                         | 2 Brush Holder                    |
| D-A2–A3.5, B2, B3 . . . 2.9 ml    | 1 Brush No. 7                     |
| 3 Enamel E-57–59 . . . . . 2.9 ml | 2 Syringe Plunger                 |
| 1 Separator. . . . . 5 ml         | 3 Instructions for use            |
| 1 Diamond Polisher. . . . . 8 g   | 1 Creation VC Workbook            |

VC Kit 2 | 820202  
Effect porcelains kit

**A colour system, a layering concept and a ceramic-like feeling:**  
The extension set offers dental technicians a wide spectrum for creating a natural look, using a variety of opalescent and translucent porcelains – even for very thin layer thicknesses. The **Flow porcelains** are a real highlight. With its unique filler technology, this fluid composite shows an extremely homogeneous distribution of ultra-fine particles within the matrix. The filler bonds perfectly to the matrix, guaranteeing a very high bonding strength. The excellent polishing qualities ensure that when this almost self-polishing material is used for long-term restorations, it will give off a shine that remains for a considerable period of time. The flow materials are outstandingly suited for bonding with the Creopal Shells, and promise a durable bond with stable colouring, thanks to the latest generation of primer (VC Composite Primer). Producing restorations in the injection technique as well as the classical layering of veneers and tabletops are only some of the possibilities for the universal use of the Flow porcelains.



- |                                   |                                   |
|-----------------------------------|-----------------------------------|
| 2 Pearl Enamel PS-1, 2 . . 2.9 ml | 2 Flow Opaque Dentine             |
| 3 Effect Enamel SI-01–03. 2.9 ml  | OD-A2, A3. . . . . 2 ml           |
| 4 Transpa                         | 2 Flow Enamel EN-L, T . . . 2 ml  |
| TI-01, 02, 04, 05 . . . . 2.9 ml  | 1 Composite Primer . . . . . 3 ml |
| 3 Neck Transpa                    | 20 Dispensing Tips                |
| HT-52–54. . . . . 2.9 ml          | 2 Syringe Plunger                 |
| 9 Flow Dentine                    | 3 Instructions for use            |
| BW, A1–3,5, B1–3, HT . . 2 ml     |                                   |

VC Kit 2 Starter | 820207  
Effect porcelains introduction kit

A colour system, a layering concept and a ceramic-like feeling with a reduced selection of the most common shades – based on kit 1. Can be completed according to the requirements. The set offers dental technicians a wide spectrum for creating a natural look, using a variety of opalescent and translucent porcelains – even for very thin layer thicknesses. The **Flow porcelains** are a real highlight.



- |                                  |                                   |
|----------------------------------|-----------------------------------|
| 3 Effect Enamel SI-01–03. 2.9 ml | 2 Flow Opaque Dentine             |
| 4 Transpa                        | OD-A2, A3. . . . . 2 ml           |
| TI-01, 02, 04, 05 . . . . 2.9 ml | 1 Flow Enamel EN-L . . . . . 2 ml |
| 3 Neck Transpa                   | 1 Composite Primer . . . . . 3 ml |
| HT-52–54. . . . . 2.9 ml         | 20 Dispensing Tips                |
| 6 Flow Dentine                   | 2 Syringe Plunger                 |
| A2–3,5, B2–3, HT . . . . . 2 ml  | 3 Instructions for use            |

VC Gingiva Kit | 820203  
Complete gingiva colours and primers

**Livelier red aesthetics:**  
The gingiva colours kit makes it easy to reproduce all areas of the gum, with a completely natural look. The porcelains are available in three viscosities (liquid, paste and gel) – for a homogeneous foundation, simple processing, mixing, and very thin layer thicknesses! The liquid opaque's greatest qualities are its extraordinary covering ability and ease of use. The paste VC gingiva materials can be used individually, mixed or in combination with other shades, while the gel modifiers are ideal for small space ratios and thin layer thicknesses, due to their viscosity.



- |                                                    |                                                         |
|----------------------------------------------------|---------------------------------------------------------|
| 4 Gingiva Opaque<br>O-G3-6 . . . . . 2.4 ml        | 1 Gingiva Translucent<br>GT light-pink . . . . . 2.4 ml |
| 4 Gingiva G3-6 . . . . . 2.9 ml                    | 1 Metal Primer . . . . . 5 ml                           |
| 6 Gingiva Modifier<br>GM-1-3, 5-7 . . . . . 2.4 ml | 1 Composite Primer . . . . . 3 ml                       |
|                                                    | 20 Dispensing Tips                                      |
|                                                    | 2 Syringe Plunger                                       |
|                                                    | 3 Instructions for use                                  |



created by Uwe Gehringer, Oral Design

VC Make up Neo Kit | 820204  
Universal stains

**Amazingly simple and varied:**  
The fine-grained, nano-filled stains and modifiers allow simple and efficient colouring, for both internal and external (transparent) characterisation of the layer structures. The high-quality VC Make up Neo porcelains are also well-suited for the palatal modification of Creopal denture teeth and shells, for greater efficiency and versatility. The special nano-filler technology helps to guarantee a perfect particle distribution, resulting in a very high resistance to abrasion and colour retention – even with thin layer thicknesses. In many cases, mechanical polishing can be omitted or clearly shortened. The two translucent gel materials, Translucent h.v. (high viscosity) and l.v. (low viscosity), are suitable for mixing with the Neo colours for internal characterisation.



- |                                                  |                                       |
|--------------------------------------------------|---------------------------------------|
| 3 Make up Neo<br>MUN-A-C . . . . . 2.6 ml        | 2 Translucent h.v., l.v. . . . 2.4 ml |
| 12 Make up Neo<br>MUN-1-4, 7-12, IL, G. . 2.6 ml | 10 Brush Round, Size 1                |
| 2 Make up Neo<br>MUN-CL, MUN-CL h.c. . 5 ml      | 10 Brush Flat, Size 1                 |
|                                                  | 2 Brush Holder                        |
|                                                  | 2 Instructions for use                |

VC Trial Kit A2 | 820208  
The easy introduction of the system

A colour system, a layering concept and a ceramic-like feeling with a reduced selection – matched to the most common shade A2. The high-quality 1-shade set offers dental technicians a wide spectrum for creating a natural look – even for very thin layer thicknesses. The **Flow porcelains** are a real highlight – also in combination with Creopal- and Creopal Shell-denture teeth.



- |                                               |                                   |
|-----------------------------------------------|-----------------------------------|
| 1 Base Opaque O-BO. . . 2.4 ml                | 1 Flow Enamel EN-L . . . . 2.9 ml |
| 1 Cervical Opaque<br>O-CO. . . . . 2.4 ml     | 1 Separator . . . . . 5 ml        |
| 1 Opaque O-A2 . . . . . 2.4 ml                | 1 Diamond Polisher . . . . . 8 g  |
| 1 Opaque Dentine<br>OD-32 . . . . . 2.9 ml    | 1 Air Protect . . . . . 10 ml     |
| 1 Dentine D-A2 . . . . . 2.9 ml               | 1 Composite Primer . . . . . 3 ml |
| 1 Enamel E-58 . . . . . 2.9 ml                | 1 Mixing Pad                      |
| 1 Make up Neo MUN-2. . 2.6 ml                 | 20 Dispensing Tips                |
| 1 Flow Dentine A2. . . . . 2 ml               | 2 Syringe Plunger                 |
| 1 Flow Opaque Dentine<br>OD-A2 . . . . . 2 ml | 4 Instructions for use            |
|                                               | 1 Creation VC Workbook            |

VC Accessory Kit | 820205  
Accessories and primers

Brushes, polishing paste or a specific primer: this kit has the right accessory for every laboratory table and indication.



- 2 Syringe Plunger

1 Separator . . . . . 5ml

1 Diamond Polisher . . . . . 8g

1 Air Protect . . . . . 10ml

1 Composite Primer . . . . . 3ml

1 Metal Primer . . . . . 5ml

1 Ceramic Primer . . . . . 3ml
- 10 Disposable Palette

10 Brush Round, Size 1

10 Brush Flat, Size 1

1 Brush No. 7

2 Brush Holder

1 Mixing Pad

4 Instructions for use

VC Single Products and Bonding Agents

The right tools for a strong bond to different materials, or a long-lasting shine.



**VC Metal Primer | 820154**  
Ensures secure bonding between the framework and the first opaque layer (VC Base Opaque).


**VC Ceramic Primer | 820156**  
This primer's special composition makes it an important part of the adhesive bond to the first opaque layer, and subsequently, the composite bond.

**VC Composite Primer | 820155**  
Light-curing primer for recreating the inhibition layer (reaction layer). Builds the bond from already polymerised material to material which has not yet been polymerised.

**VC Separator | 820151**  
Composite-plaster insulation used with the working plaster model for direct contact with the composite.

**VC Air Protect | 820153**  
Acid blocker which protects the inhibition layer after polymerisation, in order to maintain the surface structure. This guarantees complete polymerisation of the surface.

**VC Diamond Polisher | 820152**  
Extremely fine-grained diamond polishing paste which quickly and efficiently creates a controlled shine.



# DENTURE TEETH

“Without any loss of aesthetics, the natural, slightly rectangular tooth shapes can perfectly be individualised and combined also due to the skilfully arranged interdental closing ridges.

My conclusion: Creopal teeth and shells are economic and aesthetic down to the smallest detail.”

*Michael Zangl*



## CREOPAL COMPOSITE TEETH AND SHELLS – exclusive range of denture teeth: multifaceted play of colour and shape

High aesthetics, stability and efficiency: Creopal is a modern system of composite teeth with enhanced strength – the result of intercultural exchange. Thus the gnathology and the surface texture of the anterior teeth were designed in a Willi Geller project lasting several years, while the posterior teeth were developed in parallel by Yasuhiro Odanaka, among others.

The **ten anterior and five posterior tooth shapes** for the maxilla and mandibula are perfectly coordinated with each other. Particularly impressive features are their natural morphology and surface texture as well as an extraordinary gloss. As a result of the balanced filler content, the highly compacted composite teeth display high abrasion and flexural strength.

The exclusive range of artificial teeth also includes multifunctional **composite shells for the anterior and posterior area**. These Creopal Shells, despite their extremely thin layer thickness,

achieve the same colour effect and the same brightness value as full contour teeth – creating a masterful impact!

**Individually adaptable and multiple indications:** The Creopal system is available in 16 Vita® and two bleach shades, for perfect colour adaptation with acrylic and ceramic materials, the Creopal Shell system has a limited colour range. The artificial teeth produced by a four-layer process and the Creopal Shells with a two-layer method, are suitable for partial, combination and implant restorations as well as full dentures and as a set-up or mock-up for aesthetic try-ins and backward planning. The slightly rectangular tooth shapes can be perfectly individualised and combined – without any loss of aesthetics – best with Creation VC!

**Creopal teeth and shells – for individual tooth aesthetics.**

## THE BRILLIANCE OF CREOPAL AND CREOPAL SHELL:

- System-independent, lifelike denture teeth and veneer shells made from composite
- Tooth shapes with enhanced opalescence and fluorescence – for more light dynamics
- Rectangular basic forms can easily be individualised without loss of colour
- Closed interdental spaces for a natural set-up and papilla shaping
- Universal indications in temporary restorations, partial and full dentures

CREOPAL: SHELL AND FULL CONTOUR

|                 |   |   |   |   |   |  |   |   |   |   |  |   |   |   |   |  |   |   |   |   |   |
|-----------------|---|---|---|---|---|--|---|---|---|---|--|---|---|---|---|--|---|---|---|---|---|
| shell anterior  | X | X | X | X | - |  | X | X | X | - |  | - | - | - | - |  | - | - | - | - | X |
| shell posterior | X | X | X | X | - |  | X | X | X | - |  | - | - | - | - |  | - | - | - | - | - |
| full anterior   | X | X | X | X | X |  | X | X | X | X |  | X | X | X | X |  | X | X | X | X | X |
| full posterior  | X | X | X | X | X |  | X | X | X | X |  | X | X | X | X |  | X | X | X | - | - |





**Creopal Anterior**  
Fully anatomical composite anterior tooth.

**Single sets, shapes and shades**  
Creopal upper jaw anterior teeth AU, 6-piece . . . . . 1 set  
6 shapes,  
16 tooth shades A1–D4 and 2 bleach shades FW2, FW5

Creopal lower jaw anterior teeth AL, 6-piece . . . . . 1 set  
4 shapes,  
16 tooth shades A1–D4 and 2 bleach shades FW2, FW5



**Creopal Posterior**  
Fully anatomical composite posterior tooth.

**Single sets, shapes and shades**  
Creopal upper jaw posterior teeth PU, 8-piece . . . . . 1 set  
4 shapes, 16 tooth shades A1–D4

Creopal lower jaw posterior teeth PL, 8-piece . . . . . 1 set  
4 shapes, 16 tooth shades A1–D4



**Creopal Shell Anterior**  
Veneer composite anterior tooth.

**Single sets, shapes and shades**  
Creopal Shell upper jaw anterior teeth SAU 6-piece . . . . . 1 set  
4 shapes,  
7 tooth shades A1–3.5, B1–3 and 1 bleach shade FW5

Creopal Shell lower jaw anterior teeth SAL, 6-piece . . . . . 1 set  
3 shapes,  
7 tooth shades A1–3.5, B1–3 and 1 bleach shade FW5



**Creopal Shell Posterior**  
Veneer composite posterior tooth.

**Single sets, shapes and shades**  
Creopal Shell upper jaw posterior teeth SPU, 8-piece . . . . . 1 set  
4 shapes, 7 tooth shades A1–3.5, B1–3

Creopal Shell lower jaw posterior teeth SPL, 8-piece . . . . . 1 set  
4 shapes, 7 tooth shades A1–3.5, B1–3



**Creation Tooth Cabinet**  
The tooth cabinet with five drawers offers enough space for up to 300 sets (60 anterior sets or 40 posterior sets per inlay). The drawer bases are made of anodised, black moulded sheet metal parts, and can be opened by lightly touching the front. Each drawer is fitted with a acrylic glass insertion plate at the front for individual labelling. Tooth cabinet without inlays and sets. 6-piece and 8-piece inlays are available individually.  
Dimensions: 44.5 x 39.8 x 19.5 cm

Creation Tooth Cabinet | 1000300030  
Assortment inlay for 72 anterior sets | PEARL00003  
Assortment inlay for 48 posterior sets | PEARL00004



**Creopal Denture Tooth Presentation Boxes**  
High-value aluminium chests with 33 sets in all shapes and Creopal primer.

- Creopal Presentation Box . . . . . A1 | 1100000011
- Creopal Presentation Box . . . . . A2 | 1100000012
- Creopal Presentation Box . . . . . A3 | 1100000013
- Creopal Presentation Box . . . . . A3.5 | 1100000014
- Creopal Presentation Box . . . . . B2 | 1100000017
- Creopal Presentation Box . . . . . B3 | 1100000018



**Creopal Shade Guide**  
Tooth shade guide with 16 tooth shades A1–D4 and 2 bleach shades FW2 and FW5.

Creopal Shade Guide | 1100000001




**Creopal Primer**  
Bonding agent for a secure bond between prosthetic tooth and prosthetic resin.

Creopal Primer | 1200000001



Type AU – upper anterior



Type AL – lower anterior



**Contact Relationship in Statics**  
Flat intercuspation  
  
Calotte: 140  
Can readily be used in a combination technique and complete prosthetics due to the anatomical shape and the cusp inclination.





CREAPEARL 2 PMMA TEETH –  
cosmetic dentistry that works functional and looks natural

Simply because being natural matters: the most fascinating characteristic of the Creapearl 2 anterior and posterior teeth is not only their brilliant aesthetics, but also their universal function. This means that the maximum demand for individuality can be reached. Minimal grinding on the body of the tooth also makes it easier to adjust the form of the tooth.

Thanks to the elaborate manufacturing process, high-value materials and a particularly strong dentine core, this range of anterior and posterior teeth, independent of any system, guarantees dental technicians **prosthetic teeth that are 100 % true to their colour, shape and function** – adapted to the Vita® shade system.

Creapearl 2 Anterior

The lively colour effect of the Creapearl 2 Anterior makes it a classic anterior tooth with a wide variety of shapes. Its other distinguishing features are its **almost natural layering** with increased opalescence, fluorescence and translucence, and its lifelike surface structure.

Creapearl 2 Posterior

The posterior tooth has a **high level of abrasion resistance** and increased vertical stability and resilience – for a long life and resistance to plaque. Its **natural surface morphology** also makes it feel comfortable in the mouth.

The basic concept of the Creapearl 2 tooth range involves the reproduction of a natural tooth's function and character, using strengthened PMMA resin. The form follows the function, bringing the **fundamental principles of Creation** into effect:

- Raised opacity
- Chromatic dentine
- Opalescent glaze coating
- Individual high-contrast layering (coloured/white, opaque/translucent)

When the form and surface texture of the Creapearl 2 anterior teeth was conceived by master of ceramics Willi Geller, the MDT Dieter Schulz came up with the concept for the posterior teeth, using the occlusal compass and drawing on his great experience in NAT/NFR recovery technology. Thanks to this easy-to-understand positioning concept and the **multi-functional occlusal surface design** (with free spaces in the occlusion), all types of position and biting situations can be implemented.

Creapearl 2: natural and functional across all our designs!

THE BRILLIANCE OF  
CREAPEARL 2:

- Anterior and posterior tooth range with perfect colouring and function
- Dynamic colour penetration for a reduction without colour wearing away
- Closed approximal areas for a vital red/white aesthetic ("Geller Wing")
- More expressive and natural structure of the occlusal surface relief
- Outstanding resistance to wear and tear
- Strong, characteristic tooth shape – can be used by anyone, for anyone, and adapted to suit the individual

CREAPEARL 2

|           |   |   |   |   |   |  |   |   |   |   |  |   |   |   |   |  |   |   |   |
|-----------|---|---|---|---|---|--|---|---|---|---|--|---|---|---|---|--|---|---|---|
| anterior  | X | X | X | X | X |  | X | X | X | X |  | X | X | X | X |  | X | X | X |
| posterior | X | X | X | X | X |  | X | X | X | X |  | X | X | X | X |  | X | X | X |



Creation Tooth Cabinet

The tooth cabinet with five drawers offers enough space for up to 300 sets (60 anterior sets or 40 posterior sets per inlay). The drawer bases are made of anodised, black moulded sheet metal parts, and can be opened by lightly touching the front. Each drawer is fitted with a acrylic glass insertion plate at the front for individual labelling. Tooth cabinet without inlays and sets. 6-piece and 8-piece inlays are available individually. Dimensions: 44.5 x 39.8 x 19.5 cm

Creation Tooth Cabinet | 1000300030

Assortment inlay for 72 anterior sets | PEARL00003  
Assortment inlay for 48 posterior sets | PEARL00004



**Creapearl 2 Denture Tooth Presentation Boxes**  
High-value aluminium chests with 33 sets in all shapes.

Creapearl 2 Presentation Box . . . . . A1 | 1000300011  
Creapearl 2 Presentation Box . . . . . A2 | 1000300012  
Creapearl 2 Presentation Box . . . . . A3 | 1000300013  
Creapearl 2 Presentation Box . . . . . A3.5 | 1000300014



**Creapearl 2 Shade Guide**  
Tooth shade guide with 16 tooth shades A1–D4.

Creapearl 2 Shade Guide | 1000300009



**Creapearl 2 Anterior**  
Fully anatomical anterior tooth made from strengthened PMMA resin.

**Single sets, shapes and shades**  
Creapearl 2 upper jaw anterior teeth A and V, 6-piece . . . . 1 set  
10 shapes, 16 tooth shades A1–D4

Creapearl 2 lower jaw anterior teeth U, 6-piece . . . . . 1 set  
5 shapes, 16 tooth shades A1–D4



**Creapearl 2 Posterior**  
Fully anatomical posterior tooth made from strengthened PMMA resin.

**Single sets, shapes and shades**  
Creapearl 2 upper jaw posterior teeth D, E, F, G, 8-piece . . 1 set  
9 shapes, 16 tooth shades A1–D4

Creapearl 2 lower jaw posterior teeth D, E, F, G, 8-piece . . 1 set  
9 shapes, 16 tooth shades A1–D4

| Technical specifications | Dynamicline (D)                 | Balanceline (E, F)                   | Gnathosline (G)                    |
|--------------------------|---------------------------------|--------------------------------------|------------------------------------|
| Dimensions               | S, M, L                         | S, M                                 | S, M                               |
| Morphology               | Constructional morphology (NAT) | Extended functional morphology (NFR) | Natural shape and cusp inclination |


Constructional morphology (NAT), extended functional morphology (NFR), MDT Dieter Schulz



**Type A Anatomic – upper anterior teeth with natural root-shaped emergence profile**  
The denture tooth is a perfect anatomical copy of a natural tooth.



**Type V Vario – upper anterior teeth for individual contouring**  
The clinical tooth crown can be modified for an individual appearance, with or without grinding down the tooth neck, thanks to the model of the gingival edge.



**Type U – lower anterior teeth**  
Can be arranged in the same way as Type V Vario, without affecting the anatomical circumstances. In this way, the length of the tooth itself can be determined.



**Dynamicline**  
Calotte: 140  
Condyle inclination (HCN) FH 45° / CE 30°  
Bennett angle 10°, ISS 0.5mm  
Lateroretrusion (LRT) 20°  
Retrusion (RT) 0.5mm



**Balanceline E**  
Calotte: 140  
Condyle inclination (HCN) FH 30° / CE 15°  
Bennett angle 10°, ISS 1.0mm  
Lateroretrusion (LRT) 20°  
Retrusion (RT) 1.0mm



**Balanceline F**  
Basal reduced  
Calotte: 140  
Condyle inclination (HCN) FH 30° / CE 15°  
Bennett angle 10°, ISS 1.0mm  
Lateroretrusion (LRT) 20°  
Retrusion (RT) 1.0mm



**Gnathosline G / Willi Geller**  
Calotte: 140  
Due to the natural shape and the cusp inclination, the denture teeth can be used for partial, combined and full denture restorations.



WILLI GELLER  
*Creation*  
SIMPLY BRILLIANT

Creapearl 2

WILLI GELLER  
*Creation*  
SIMPLY BRILLIANT

CONTENT

CREATION CERAMICS AND ACCESSORIES

All kits at a glance ..... 3

CC - Classic Metal Ceramic ..... 4

ZI-CT - Zirconium Dioxide Ceramic..... 13

ZI-F - Low-Fusing Zirconium Dioxide Ceramic ..... 17

LS/LS Press - Lithium Disilicate  
Veneering and Pressable Ceramic ..... 21

CreaColor – Stains und Modifiers ..... 23

Creation Liquids ..... 24

Ceramic Accessories ..... 24

Creation Shade Guides ..... 25

Aqualine ..... 26

CREATION COMPOSITE

VC - Light-Curing Composite ..... 27

CREATION DENTURE TEETH

Creopal ..... 31

Creapearl 2 ..... 41

Printing errors excepted. Subject to change at any time. The general terms and conditions of Creation Willi Geller International GmbH apply (www.creation-willigeller.com)

Shades A1–D4 are referring to the Vita® Shade Guide. Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

| Article no. | Article name                                             | Content |
|-------------|----------------------------------------------------------|---------|
|             | CERAMIC                                                  |         |
| 000300      | CC / Kit Creapast                                        | 1 Kit   |
| 000800      | CC / Powder Opaque Kit (inkl. Shade Guide 000068)        | 1 Kit   |
| 000804      | CC / Kit 1 804 (inkl. Shade Guide 000064)                | 1 Kit   |
| 001508      | CC / Kit 2 1508 (inkl. Shade Guide 000073)               | 1 Kit   |
| 000915      | CC / Starter Kit                                         | 1 Kit   |
| 000920      | CC / Introduction Kit                                    | 1 Kit   |
| 001590      | CC / Gingiva Kit (inkl. Shade Guide 000048)              | 1 Kit   |
| 000250      | CC / Bleach Shades Kit (inkl. Shade Guide 000049)        | 1 Kit   |
| 001503      | CC / Shoulder Kit (inkl. Shade Guide 000047)             | 1 Kit   |
| 000801      | CC / Kit 1 Transition Dentine (inkl. Shade Guide 000087) | 1 Kit   |
| 620501      | ZI-CT / Kit 1                                            | 1 Kit   |
| 620502      | ZI-CT / Kit 2                                            | 1 Kit   |
| 620503      | ZI-CT / Starter Kit                                      | 1 Kit   |
| 620506      | ZI-CT / Introduction Kit                                 | 1 Kit   |
| 620504      | ZI-CT / Gingiva Kit                                      | 1 Kit   |
| 620505      | ZI-CT / Bleach Kit                                       | 1 Kit   |
| 620507      | ZI-CT / Shoulder Kit                                     | 1 Kit   |
| 653010      | ZI-F / Kit 1 (inkl. Shade Guide 000076)                  | 1 Kit   |
| 653015      | ZI-F / Kit 2 (inkl. Shade Guide 000077)                  | 1 Kit   |
| 653025      | ZI-F / Starter Kit                                       | 1 Kit   |
| 653030      | ZI-F / Introduction Kit                                  | 1 Kit   |
| 603030      | ZI-F / Gingiva Kit (inkl. Shade Guide 000048)            | 1 Kit   |
| 652050      | ZI-F / Bleach Shades Kit (inkl. Shade Guide 000049)      | 1 Kit   |
| 830080      | LS Press / Starter Kit                                   | 1 Kit   |
| 830081      | LS Press / Introduction Kit A1, A2                       | 1 Kit   |
| 830082      | LS Press / Introduction Kit A3                           | 1 Kit   |
| 830083      | LS / Kit                                                 | 1 Kit   |
| 001200      | CreaColor Kit                                            | 1 Kit   |
| 001100      | Make Up Neo Kit                                          | 1 Kit   |
| 001150      | In Nova Neo Kit                                          | 1 Kit   |
|             | COMPOSITE                                                |         |
| 820201      | VC / Kit 1                                               | 1 Kit   |
| 820202      | VC / Kit 2                                               | 1 Kit   |
| 820203      | VC / Gingiva Kit                                         | 1 Kit   |
| 820204      | VC / Make up Neo Kit                                     | 1 Kit   |
| 820205      | VC / Accessory-Kit                                       | 1 Kit   |
| 820206      | VC / Kit 1 Starter                                       | 1 Kit   |
| 820207      | VC / Kit 2 Starter                                       | 1 Kit   |
| 820208      | VC / Trial Kit                                           | 1 Kit   |

CREATION CC | CLASSIC METAL CERAMIC

| Article no. | Article name               | Content |
|-------------|----------------------------|---------|
|             | <b>PASTES 3 g – 18 g</b>   |         |
| 000080      | CC.LF/Crea Alloy Bond      | 2 x 4 g |
| 000301      | CC/Creapast A1             | 3 g     |
| 000302      | CC/Creapast A2             | 3 g     |
| 000303      | CC/Creapast A3             | 3 g     |
| 000304      | CC/Creapast A3,5           | 3 g     |
| 000305      | CC/Creapast A4             | 3 g     |
| 000306      | CC/Creapast B1             | 3 g     |
| 000307      | CC/Creapast B2             | 3 g     |
| 000308      | CC/Creapast B3             | 3 g     |
| 000309      | CC/Creapast B4             | 3 g     |
| 000310      | CC/Creapast C1             | 3 g     |
| 000311      | CC/Creapast C2             | 3 g     |
| 000312      | CC/Creapast C3             | 3 g     |
| 000313      | CC/Creapast C4             | 3 g     |
| 000314      | CC/Creapast D2             | 3 g     |
| 000315      | CC/Creapast D3             | 3 g     |
| 000316      | CC/Creapast D4             | 3 g     |
| 000317      | CC/Creapast IO white       | 3 g     |
| 000318      | CC/Creapast IO violet      | 3 g     |
| 000319      | CC/Creapast IO ochre       | 3 g     |
| 000320      | CC/Creapast IO gold yellow | 3 g     |
| 000322      | CC/Creapast WOP sand       | 3 g     |
| 000323      | CC/Creapast WOP sahara     | 3 g     |
| 000324      | CC/Creapast WOP gold       | 3 g     |
| 000325      | CC/Creapast G gingiva      | 3 g     |
| 000326      | CC/Bleach Creapast O-AB    | 3 g     |
| 006301      | CC/Creapast A1             | 18 g    |
| 006302      | CC/Creapast A2             | 18 g    |
| 006303      | CC/Creapast A3             | 18 g    |
| 006304      | CC/Creapast A3,5           | 18 g    |
| 006305      | CC/Creapast A4             | 18 g    |
| 006306      | CC/Creapast B1             | 18 g    |
| 006307      | CC/Creapast B2             | 18 g    |
| 006308      | CC/Creapast B3             | 18 g    |
| 006309      | CC/Creapast B4             | 18 g    |
| 006310      | CC/Creapast C1             | 18 g    |
| 006311      | CC/Creapast C2             | 18 g    |
| 006312      | CC/Creapast C3             | 18 g    |
| 006313      | CC/Creapast C4             | 18 g    |
| 006314      | CC/Creapast D2             | 18 g    |
| 006315      | CC/Creapast D3             | 18 g    |
| 006316      | CC/Creapast D4             | 18 g    |
| 006317      | CC/Creapast IO white       | 18 g    |
| 006318      | CC/Creapast IO violet      | 18 g    |
| 006319      | CC/Creapast IO ochre       | 18 g    |
| 006320      | CC/Creapast IO gold yellow | 18 g    |
| 006322      | CC/Creapast WOP sand       | 18 g    |

CREATION CC | CLASSIC METAL CERAMIC

| Article no. | Article name                            | Content |
|-------------|-----------------------------------------|---------|
| 006323      | CC/Creapast WOP sahara                  | 18 g    |
| 006324      | CC/Creapast WOP gold                    | 18 g    |
| 006325      | CC/Creapast G gingiva                   | 18 g    |
| 006326      | CC/Bleach Creapast O-AB                 | 18 g    |
|             | <b>PORCELAIN 10 g – 15 g</b>            |         |
| 000174      | CC/Glaze GL                             | 10 g    |
| 001511      | CC/Opaque Modifier OM-11 white          | 15 g    |
| 001512      | CC/Opaque Modifier OM-12 pink           | 15 g    |
| 001513      | CC/Opaque Modifier OM-13 lilac          | 15 g    |
| 001514      | CC/Opaque Modifier OM-14 grey           | 15 g    |
| 001515      | CC/Opaque Modifier OM-15 sahara         | 15 g    |
| 001516      | CC/Opaque Modifier OM-16 ocker          | 15 g    |
| 001517      | CC/Opaque Modifier OM-17 brown          | 15 g    |
| 001518      | CC/Opaque Modifier OM-G gingiva         | 15 g    |
| 001521      | CC/Shoulder SP-21 neutral               | 15 g    |
| 001522      | CC/Shoulder SP-22 flamingo              | 15 g    |
| 001523      | CC/Shoulder SP-23 sand                  | 15 g    |
| 001524      | CC/Shoulder SP-24 gold                  | 15 g    |
| 001525      | CC/Shoulder SP-25 red brown             | 15 g    |
| 001526      | CC/Shoulder SP-26 yellow grey           | 15 g    |
| 001527      | CC/Shoulder SP-27 ivory                 | 15 g    |
| 001528      | CC/Shoulder SP-28 opaque yellow         | 15 g    |
| 001529      | CC/Shoulder SP-29 honey                 | 15 g    |
| 001530      | CC/Shoulder SP-G gingiva                | 15 g    |
| 000256      | CC/Bleach Shoulder SP-AB                | 15 g    |
| 001532      | CC/Opaque Dentine OD-32 (A2-A3)         | 15 g    |
| 001537      | CC/Opaque Dentine OD-37 (B4)            | 15 g    |
| 001541      | CC/Opaque Dentine OD-41 (D3)            | 15 g    |
| 001543      | CC/Opaque Dentine OD-43 (ivory)         | 15 g    |
| 001544      | CC/Opaque Dentine OD-44 (cuba)          | 15 g    |
| 001551      | CC/Neck Transpa HT-51 iridescent        | 15 g    |
| 001552      | CC/Neck Transpa HT-52 khaki             | 15 g    |
| 001553      | CC/Neck Transpa HT-53 sand              | 15 g    |
| 001554      | CC/Neck Transpa HT-54 honey             | 15 g    |
| 001555      | CC/Neck Transpa HT-55 bordeaux          | 15 g    |
| 001556      | CC/Neck Transpa HT-56 ochre             | 15 g    |
| 001561      | CC/Make In MI-61 ivory                  | 15 g    |
| 001562      | CC/Make In MI-62 lemon                  | 15 g    |
| 001563      | CC/Make In MI-63 honey                  | 15 g    |
| 001564      | CC/Make In MI-64 flamingo               | 15 g    |
| 001565      | CC/Make In MI-65 gold                   | 15 g    |
| 001566      | CC/Make In MI-66 olive                  | 15 g    |
| 001568      | CC/Approximal Dentine AD-1 light yellow | 15 g    |
| 001569      | CC/Approximal Dentine AD-2 orange       | 15 g    |
| 001570      | CC/Opaque Dentine OD-A1                 | 15 g    |
| 001571      | CC/Opaque Dentine OD-A2                 | 15 g    |
| 001572      | CC/Opaque Dentine OD-A3                 | 15 g    |

CREATION CC | CLASSIC METAL CERAMIC

| Article no. | Article name                   | Content |
|-------------|--------------------------------|---------|
| 001573      | CC/Opaque Dentine OD-A3,5      | 15 g    |
| 001574      | CC/Opaque Dentine OD-A4        | 15 g    |
| 001575      | CC/Opaque Dentine OD-B1        | 15 g    |
| 001576      | CC/Opaque Dentine OD-B2        | 15 g    |
| 001577      | CC/Opaque Dentine OD-B3        | 15 g    |
| 001578      | CC/Opaque Dentine OD-B4        | 15 g    |
| 001579      | CC/Opaque Dentine OD-C1        | 15 g    |
| 001580      | CC/Opaque Dentine OD-C2        | 15 g    |
| 001581      | CC/Opaque Dentine OD-C3        | 15 g    |
| 001582      | CC/Opaque Dentine OD-C4        | 15 g    |
| 001583      | CC/Opaque Dentine OD-D2        | 15 g    |
| 001584      | CC/Opaque Dentine OD-D3        | 15 g    |
| 001585      | CC/Opaque Dentine OD-D4        | 15 g    |
| 001591      | CC/Gingiva G1 purple           | 15 g    |
| 001592      | CC/Gingiva G2 dark pink        | 15 g    |
| 001593      | CC/Gingiva G3 light pink       | 15 g    |
| 001594      | CC/Gingiva G4 flamingo         | 15 g    |
| 001595      | CC/Gingiva G5 rose             | 15 g    |
| 001596      | CC/Gingiva G6 dark pink opaque | 15 g    |
| 001597      | CC/Gingiva GN neutral          | 15 g    |
| 001598      | CC/Gingiva G7 violet           | 15 g    |
|             | <b>PORCELAIN 20 g</b>          |         |
| 000201      | CC/Opaque O-A1                 | 20 g    |
| 000202      | CC/Opaque O-A2                 | 20 g    |
| 000203      | CC/Opaque O-A3                 | 20 g    |
| 000204      | CC/Opaque O-A3,5               | 20 g    |
| 000205      | CC/Opaque O-A4                 | 20 g    |
| 000206      | CC/Opaque O-B1                 | 20 g    |
| 000207      | CC/Opaque O-B2                 | 20 g    |
| 000208      | CC/Opaque O-B3                 | 20 g    |
| 000209      | CC/Opaque O-B4                 | 20 g    |
| 000210      | CC/Opaque O-C1                 | 20 g    |
| 000211      | CC/Opaque O-C2                 | 20 g    |
| 000212      | CC/Opaque O-C3                 | 20 g    |
| 000213      | CC/Opaque O-C4                 | 20 g    |
| 000214      | CC/Opaque O-D2                 | 20 g    |
| 000215      | CC/Opaque O-D3                 | 20 g    |
| 000216      | CC/Opaque O-D4                 | 20 g    |
| 000217      | CC/Opaque FOP                  | 20 g    |
| 000231      | CC/Dentine D-A1                | 20 g    |
| 000232      | CC/Dentine D-A2                | 20 g    |
| 000233      | CC/Dentine D-A3                | 20 g    |
| 000234      | CC/Dentine D-A3,5              | 20 g    |
| 000235      | CC/Dentine D-A4                | 20 g    |
| 000236      | CC/Dentine D-B1                | 20 g    |
| 000237      | CC/Dentine D-B2                | 20 g    |
| 000238      | CC/Dentine D-B3                | 20 g    |

CREATION CC | CLASSIC METAL CERAMIC

| Article no. | Article name                         | Content |
|-------------|--------------------------------------|---------|
| 000239      | CC/Dentine D-B4                      | 20 g    |
| 000240      | CC/Dentine D-C1                      | 20 g    |
| 000241      | CC/Dentine D-C2                      | 20 g    |
| 000242      | CC/Dentine D-C3                      | 20 g    |
| 000243      | CC/Dentine D-C4                      | 20 g    |
| 000244      | CC/Dentine D-D2                      | 20 g    |
| 000245      | CC/Dentine D-D3                      | 20 g    |
| 000246      | CC/Dentine D-D4                      | 20 g    |
| 000247      | CC/Dentine D-CW California white     | 20 g    |
| 000251      | CC/Bleach Opaque O-AB                | 20 g    |
| 000252      | CC/Bleach Dentine BD-A               | 20 g    |
| 000253      | CC/Bleach Dentine BD-B               | 20 g    |
| 000254      | CC/Bleach Dentine BD-B0              | 20 g    |
| 000255      | CC/Bleach Enamel S-AB                | 20 g    |
| 000261      | CC/Enamel S-057                      | 20 g    |
| 000262      | CC/Enamel S-058                      | 20 g    |
| 000263      | CC/Enamel S-059                      | 20 g    |
| 000264      | CC/Enamel S-060                      | 20 g    |
| 000265      | CC/Effect Enamel SI-01 light yellow  | 20 g    |
| 000266      | CC/Effect Enamel SI-02 medium yellow | 20 g    |
| 000267      | CC/Effect Enamel SI-03 heavy yellow  | 20 g    |
| 000268      | CC/Effect Enamel SI-04 light orange  | 20 g    |
| 000269      | CC/Effect Enamel SI-05 medium orange | 20 g    |
| 000270      | CC/Effect Enamel SI-06 heavy orange  | 20 g    |
| 000271      | CC/Window UC                         | 20 g    |
| 000272      | CC/Transpa NT neutral                | 20 g    |
| 000273      | CC/Clear CL-O                        | 20 g    |
| 000274      | CC/Glaze GL                          | 20 g    |
| 000275      | CC/Transpa OT opal                   | 20 g    |
| 000276      | CC/Transpa TI-1 blue                 | 20 g    |
| 000277      | CC/Transpa TI-2 white                | 20 g    |
| 000278      | CC/Transpa TI-3 pink                 | 20 g    |
| 000279      | CC/Transpa TI-4 yellow               | 20 g    |
| 000280      | CC/Transpa TI-5 grey                 | 20 g    |
| 000281      | CC/Pearl Enamel PS-0 white           | 20 g    |
| 000282      | CC/Pearl Enamel PS-1 blue/white      | 20 g    |
| 000283      | CC/Pearl Enamel PS-2 flamingo        | 20 g    |
| 000284      | CC/Pearl Enamel PS-3 orange          | 20 g    |
| 000285      | CC/Effect Enamel SO-10 blue          | 20 g    |
| 000286      | CC/Effect Enamel SO-11 orange        | 20 g    |
| 000290      | CC/Correction Powder KM              | 20 g    |
| 003201      | CC/Transition Dentine TD-A1          | 20 g    |
| 003202      | CC/Transition Dentine TD-A2          | 20 g    |
| 003203      | CC/Transition Dentine TD-A3          | 20 g    |
| 003204      | CC/Transition Dentine TD-A3,5        | 20 g    |
| 003205      | CC/Transition Dentine TD-A4          | 20 g    |
| 003206      | CC/Transition Dentine TD-B1          | 20 g    |
| 003207      | CC/Transition Dentine TD-B2          | 20 g    |

CREATION CC | CLASSIC METAL CERAMIC

| Article no. | Article name                     | Content |
|-------------|----------------------------------|---------|
| 003208      | CC/Transition Dentine TD-B3      | 20 g    |
| 003209      | CC/Transition Dentine TD-B4      | 20 g    |
| 003210      | CC/Transition Dentine TD-C1      | 20 g    |
| 003211      | CC/Transition Dentine TD-C2      | 20 g    |
| 003212      | CC/Transition Dentine TD-C3      | 20 g    |
| 003213      | CC/Transition Dentine TD-C4      | 20 g    |
| 003214      | CC/Transition Dentine TD-D2      | 20 g    |
| 003215      | CC/Transition Dentine TD-D3      | 20 g    |
| 003216      | CC/Transition Dentine TD-D4      | 20 g    |
| 003217      | CC/Transition Dentine TD-BA      | 20 g    |
|             |                                  |         |
|             | <b>PORCELAIN 50 g</b>            |         |
| 000501      | CC/Opaque O-A1                   | 50 g    |
| 000502      | CC/Opaque O-A2                   | 50 g    |
| 000503      | CC/Opaque O-A3                   | 50 g    |
| 000504      | CC/Opaque O-A3,5                 | 50 g    |
| 000505      | CC/Opaque O-A4                   | 50 g    |
| 000506      | CC/Opaque O-B1                   | 50 g    |
| 000507      | CC/Opaque O-B2                   | 50 g    |
| 000508      | CC/Opaque O-B3                   | 50 g    |
| 000509      | CC/Opaque O-B4                   | 50 g    |
| 000510      | CC/Opaque O-C1                   | 50 g    |
| 000511      | CC/Opaque O-C2                   | 50 g    |
| 000512      | CC/Opaque O-C3                   | 50 g    |
| 000513      | CC/Opaque O-C4                   | 50 g    |
| 000514      | CC/Opaque O-D2                   | 50 g    |
| 000515      | CC/Opaque O-D3                   | 50 g    |
| 000516      | CC/Opaque O-D4                   | 50 g    |
| 000517      | CC/Opaque O-FOP                  | 50 g    |
| 000531      | CC/Dentine D-A1                  | 50 g    |
| 000532      | CC/Dentine D-A2                  | 50 g    |
| 000533      | CC/Dentine D-A3                  | 50 g    |
| 000534      | CC/Dentine D-A3,5                | 50 g    |
| 000535      | CC/Dentine D-A4                  | 50 g    |
| 000536      | CC/Dentine D-B1                  | 50 g    |
| 000537      | CC/Dentine D-B2                  | 50 g    |
| 000538      | CC/Dentine D-B3                  | 50 g    |
| 000539      | CC/Dentine D-B4                  | 50 g    |
| 000540      | CC/Dentine D-C1                  | 50 g    |
| 000541      | CC/Dentine D-C2                  | 50 g    |
| 000542      | CC/Dentine D-C3                  | 50 g    |
| 000543      | CC/Dentine D-C4                  | 50 g    |
| 000544      | CC/Dentine D-D2                  | 50 g    |
| 000545      | CC/Dentine D-D3                  | 50 g    |
| 000546      | CC/Dentine D-D4                  | 50 g    |
| 000547      | CC/Dentine D-CW California white | 50 g    |
| 000551      | CC/Bleach Opaque O-AB            | 50 g    |
| 000552      | CC/Bleach Dentine BD-A           | 50 g    |

CREATION CC | CLASSIC METAL CERAMIC

| Article no. | Article name                         | Content |
|-------------|--------------------------------------|---------|
| 000553      | CC/Bleach Dentine BD-B               | 50 g    |
| 000554      | CC/Bleach Dentine BD-B0              | 50 g    |
| 000555      | CC/Bleach Enamel S-AB                | 50 g    |
| 000556      | CC/Bleach Shoulder SP-AB             | 50 g    |
| 000561      | CC/Enamel S-057                      | 50 g    |
| 000562      | CC/Enamel S-058                      | 50 g    |
| 000563      | CC/Enamel S-059                      | 50 g    |
| 000564      | CC/Enamel S-060                      | 50 g    |
| 000565      | CC/Effect Enamel SI-01 light yellow  | 50 g    |
| 000566      | CC/Effect Enamel SI-02 medium yellow | 50 g    |
| 000567      | CC/Effect Enamel SI-03 heavy yellow  | 50 g    |
| 000568      | CC/Effect Enamel SI-04 light orange  | 50 g    |
| 000569      | CC/Effect Enamel SI-05 medium orange | 50 g    |
| 000570      | CC/Effect Enamel SI-06 heavy orange  | 50 g    |
| 000571      | CC/Window UC                         | 50 g    |
| 000572      | CC/Transpa NT neutral                | 50 g    |
| 000573      | CC/Clear CL-O                        | 50 g    |
| 000575      | CC/Transpa OT opal                   | 50 g    |
| 000576      | CC/Transpa TI-1 blue                 | 50 g    |
| 000577      | CC/Transpa TI-2 white                | 50 g    |
| 000578      | CC/Transpa TI-3 pink                 | 50 g    |
| 000579      | CC/Transpa TI-4 yellow               | 50 g    |
| 000580      | CC/Transpa TI-5 grey                 | 50 g    |
| 000581      | CC/Pearl Enamel PS-0 white           | 50 g    |
| 000582      | CC/Pearl Enamel PS-1 blue/white      | 50 g    |
| 000583      | CC/Pearl Enamel PS-2 flamingo        | 50 g    |
| 000584      | CC/Pearl Enamel PS-3 orange          | 50 g    |
| 000585      | CC/Effect Enamel SO-10 blue          | 50 g    |
| 000586      | CC/Effect Enamel SO-11 orange        | 50 g    |
| 003501      | CC/Transition Dentine TD-A1          | 50 g    |
| 003502      | CC/Transition Dentine TD-A2          | 50 g    |
| 003503      | CC/Transition Dentine TD-A3          | 50 g    |
| 003504      | CC/Transition Dentine TD-A3,5        | 50 g    |
| 003505      | CC/Transition Dentine TD-A4          | 50 g    |
| 003506      | CC/Transition Dentine TD-B1          | 50 g    |
| 003507      | CC/Transition Dentine TD-B2          | 50 g    |
| 003508      | CC/Transition Dentine TD-B3          | 50 g    |
| 003509      | CC/Transition Dentine TD-B4          | 50 g    |
| 003510      | CC/Transition Dentine TD-C1          | 50 g    |
| 003511      | CC/Transition Dentine TD-C2          | 50 g    |
| 003512      | CC/Transition Dentine TD-C3          | 50 g    |
| 003513      | CC/Transition Dentine TD-C4          | 50 g    |
| 003514      | CC/Transition Dentine TD-D2          | 50 g    |
| 003515      | CC/Transition Dentine TD-D3          | 50 g    |
| 003516      | CC/Transition Dentine TD-D4          | 50 g    |
| 003517      | CC/Transition Dentine TD-BA          | 50 g    |
| 005011      | CC/Opaque Modifier OM-11 white       | 50 g    |
| 005012      | CC/Opaque Modifier OM-12 pink        | 50 g    |

CREATION CC | CLASSIC METAL CERAMIC

| Article no. | Article name                            | Content |
|-------------|-----------------------------------------|---------|
| 005013      | CC/Opaque Modifier OM-13 lilac          | 50 g    |
| 005014      | CC/Opaque Modifier OM-14 grey           | 50 g    |
| 005015      | CC/Opaque Modifier OM-15 sahara         | 50 g    |
| 005016      | CC/Opaque Modifier OM-16 ocker          | 50 g    |
| 005017      | CC/Opaque Modifier OM-17 brown          | 50 g    |
| 005018      | CC/Opaque Modifier OM-G gingiva         | 50 g    |
| 005021      | CC/Shoulder SP-21 neutral               | 50 g    |
| 005022      | CC/Shoulder SP-22 flamingo              | 50 g    |
| 005023      | CC/Shoulder SP-23 sand                  | 50 g    |
| 005024      | CC/Shoulder SP-24 gold                  | 50 g    |
| 005025      | CC/Shoulder SP-25 red brown             | 50 g    |
| 005026      | CC/Shoulder SP-26 yellow grey           | 50 g    |
| 005027      | CC/Shoulder SP-27 ivory                 | 50 g    |
| 005028      | CC/Shoulder SP-28 opaque yellow         | 50 g    |
| 005029      | CC/Shoulder SP-29 honey                 | 50 g    |
| 005030      | CC/Shoulder SP-G gingiva                | 50 g    |
| 005032      | CC /Opaque Dentine OD-32 (A2-A3)        | 50 g    |
| 005037      | CC /Opaque Dentine OD-37 (B4)           | 50 g    |
| 005041      | CC /Opaque Dentine OD-41 (D3)           | 50 g    |
| 005043      | CC /Opaque Dentine OD-43 (ivory)        | 50 g    |
| 005044      | CC /Opaque Dentine OD-44 (cuba)         | 50 g    |
| 005051      | CC /Neck Transpa HT-51 iridescent       | 50 g    |
| 005052      | CC/Neck Transpa HT-52 khaki             | 50 g    |
| 005053      | CC/Neck Transpa HT-53 sand              | 50 g    |
| 005054      | CC/Neck Transpa HT-54 honey             | 50 g    |
| 005055      | CC/Neck Transpa HT-55 bordeaux          | 50 g    |
| 005056      | CC/Neck Transpa HT-56 ochre             | 50 g    |
| 005061      | CC/Make In MI-61 ivory                  | 50 g    |
| 005062      | CC/Make In MI-62 lemon                  | 50 g    |
| 005063      | CC/Make In MI-63 honey                  | 50 g    |
| 005064      | CC/Make In MI-64 flamingo               | 50 g    |
| 005065      | CC/Make In MI-65 gold                   | 50 g    |
| 005066      | CC/Make In MI-66 olive                  | 50 g    |
| 005068      | CC/Approximal Dentine AD-1 light yellow | 50 g    |
| 005069      | CC/Approximal Dentine AD-2 orange       | 50 g    |
| 005070      | CC/Opaque Dentine OD-A1                 | 50 g    |
| 005071      | CC/Opaque Dentine OD-A2                 | 50 g    |
| 005072      | CC/Opaque Dentine OD-A3                 | 50 g    |
| 005073      | CC/Opaque Dentine OD-A3,5               | 50 g    |
| 005074      | CC/Opaque Dentine OD-A4                 | 50 g    |
| 005075      | CC/Opaque Dentine OD-B1                 | 50 g    |
| 005076      | CC/Opaque Dentine OD-B2                 | 50 g    |
| 005077      | CC/Opaque Dentine OD-B3                 | 50 g    |
| 005078      | CC/Opaque Dentine OD-B4                 | 50 g    |
| 005079      | CC/Opaque Dentine OD-C1                 | 50 g    |
| 005080      | CC/Opaque Dentine OD-C2                 | 50 g    |
| 005081      | CC/Opaque Dentine OD-C3                 | 50 g    |
| 005082      | CC/Opaque Dentine OD-C4                 | 50 g    |

CREATION CC | CLASSIC METAL CERAMIC

| Article no. | Article name                   | Content |
|-------------|--------------------------------|---------|
| 005083      | CC/Opaque Dentine OD-D2        | 50 g    |
| 005084      | CC/Opaque Dentine OD-D3        | 50 g    |
| 005085      | CC/Opaque Dentine OD-D4        | 50 g    |
| 005591      | CC/Gingiva G1 purple           | 50 g    |
| 005592      | CC/Gingiva G2 dark pink        | 50 g    |
| 005593      | CC/Gingiva G3 light pink       | 50 g    |
| 005594      | CC/Gingiva G4 flamingo         | 50 g    |
| 005595      | CC/Gingiva G5 rose             | 50 g    |
| 005596      | CC/Gingiva G6 dark pink opaque | 50 g    |
| 005597      | CC/Gingiva GN neutral          | 50 g    |
| 005598      | CC/Gingiva G7 violet           | 50 g    |
|             | <b>PORCELAIN 250 g</b>         |         |
| 002501      | CC/Opaque O-A1                 | 250 g   |
| 002502      | CC/Opaque O-A2                 | 250 g   |
| 002503      | CC/Opaque O-A3                 | 250 g   |
| 002504      | CC/Opaque O-A3,5               | 250 g   |
| 002505      | CC/Opaque O-A4                 | 250 g   |
| 002506      | CC/Opaque O-B1                 | 250 g   |
| 002507      | CC/Opaque O-B2                 | 250 g   |
| 002508      | CC/Opaque O-B3                 | 250 g   |
| 002509      | CC/Opaque O-B4                 | 250 g   |
| 002510      | CC/Opaque O-C1                 | 250 g   |
| 002511      | CC/Opaque O-C2                 | 250 g   |
| 002512      | CC/Opaque O-C3                 | 250 g   |
| 002513      | CC/Opaque O-C4                 | 250 g   |
| 002514      | CC/Opaque O-D2                 | 250 g   |
| 002515      | CC/Opaque O-D3                 | 250 g   |
| 002516      | CC/Opaque O-D4                 | 250 g   |
| 002531      | CC/Dentine D-A1                | 250 g   |
| 002532      | CC/Dentine D-A2                | 250 g   |
| 002533      | CC/Dentine D-A3                | 250 g   |
| 002534      | CC/Dentine D-A3,5              | 250 g   |
| 002535      | CC/Dentine D-A4                | 250 g   |
| 002536      | CC/Dentine D-B1                | 250 g   |
| 002537      | CC/Dentine D-B2                | 250 g   |
| 002538      | CC/Dentine D-B3                | 250 g   |
| 002539      | CC/Dentine D-B4                | 250 g   |
| 002540      | CC/Dentine D-C1                | 250 g   |
| 002541      | CC/Dentine D-C2                | 250 g   |
| 002542      | CC/Dentine D-C3                | 250 g   |
| 002543      | CC/Dentine D-C4                | 250 g   |
| 002544      | CC/Dentine D-D2                | 250 g   |
| 002545      | CC/Dentine D-D3                | 250 g   |
| 002546      | CC/Dentine D-D4                | 250 g   |
| 002547      | CC/California white            | 250 g   |
| 002561      | CC/Enamel S-057                | 250 g   |
| 002562      | CC/Enamel S-058                | 250 g   |

CREATION CC | CLASSIC METAL CERAMIC

| Article no. | Article name                         | Content |
|-------------|--------------------------------------|---------|
| 002563      | CC/Enamel S-059                      | 250 g   |
| 002564      | CC/Enamel S-060                      | 250 g   |
| 002565      | CC/Effect Enamel SI-01 light yellow  | 250 g   |
| 002566      | CC/Effect Enamel SI-02 medium yellow | 250 g   |
| 002567      | CC/Effect Enamel SI-03 heavy yellow  | 250 g   |
| 002568      | CC/Effect Enamel SI-04 light orange  | 250 g   |
| 002569      | CC/Effect Enamel SI-05 medium orange | 250 g   |
| 002570      | CC/Effect Enamel SI-06 heavy orange  | 250 g   |
| 002572      | CC/Transpa NT neutral                | 250 g   |
| 002573      | CC/Clear CL-O                        | 250 g   |
| 002575      | CC/Transpa OT opal                   | 250 g   |
| 002576      | CC/Transpa TI-1 blue                 | 250 g   |
| 002577      | CC/Transpa TI-2 white                | 250 g   |
| 002578      | CC/Transpa TI-3 pink                 | 250 g   |
| 002579      | CC/Transpa TI-4 yellow               | 250 g   |
| 002580      | CC/Transpa TI-5 grey                 | 250 g   |
| 002581      | CC/Pearl Enamel PS-0 white           | 250 g   |
| 002582      | CC/Pearl Enamel PS-1 blue/white      | 250 g   |
| 002583      | CC/Pearl Enamel PS-2 flamingo        | 250 g   |
| 002584      | CC/Pearl Enamel PS-3 orange          | 250 g   |
| 025071      | CC/Opaque Dentine OD-A2              | 250 g   |
| 025072      | CC/Opaque Dentine OD-A3              | 250 g   |
| 025073      | CC/Opaque Dentine OD-A3,5            | 250 g   |
| 025076      | CC/Opaque Dentine OD-B2              | 250 g   |
| 025077      | CC/Opaque Dentine OD-B3              | 250 g   |
|             | LIQUIDS                              |         |
| 000012      | Modelling Liquid high-fusing ML-H    | 25 ml   |
| 000013      | Modelling Liquid high-fusing ML-H    | 50 ml   |
| 000014      | Modelling Liquid high-fusing ML-H    | 250 ml  |
| 000015      | Modelling Liquid high-fusing ML-H    | 1000 ml |
| 000019      | Opaque Liquid high-fusing OL-H       | 25 ml   |
| 000020      | Opaque Liquid high-fusing OL-H       | 50 ml   |
| 000021      | Opaque Liquid high-fusing OL-H       | 250 ml  |
| 000010      | Universal & Glaze Liquid UF/GL       | 25 ml   |
| 000011      | Universal & Glaze Liquid UF/GL       | 50 ml   |
| 000025      | Shoulder Liquid SPL                  | 25 ml   |
| 000026      | Shoulder Liquid SPL                  | 50 ml   |
| 000090      | Moisturing Fluid                     | 25 ml   |

CREATION ZI-CT | ZIRCONIUM DIOXIDE CERAMIC

| Article no. | Article name                               | Content |
|-------------|--------------------------------------------|---------|
|             | PORCELAIN 4 g – 15 g                       |         |
| 620061      | ZI-CT/Frame Shade NT                       | 4 g     |
| 620074      | ZI-CT/Glaze GL                             | 10 g    |
| 620101      | ZI-CT/Gingiva G1 purple                    | 15 g    |
| 620102      | ZI-CT/Gingiva G2 dark pink                 | 15 g    |
| 620103      | ZI-CT/Gingiva G3 light pink                | 15 g    |
| 620104      | ZI-CT/Gingiva G4 flamingo                  | 15 g    |
| 620105      | ZI-CT/Gingiva G5 rose                      | 15 g    |
| 620106      | ZI-CT/Gingiva G6 dark pink opaque          | 15 g    |
| 620107      | ZI-CT/Gingiva G7 violet                    | 15 g    |
| 620108      | ZI-CT/Gingiva GN neutral                   | 15 g    |
| 620120      | ZI-CT/Shoulder SP-AB                       | 15 g    |
| 620121      | ZI-CT/Shoulder SP-21 neutral               | 15 g    |
| 620122      | ZI-CT/Shoulder SP-22 flamingo              | 15 g    |
| 620123      | ZI-CT/Shoulder SP-23 sand                  | 15 g    |
| 620124      | ZI-CT/Shoulder SP-24 gold                  | 15 g    |
| 620125      | ZI-CT/Shoulder SP-25 red brown             | 15 g    |
| 620126      | ZI-CT/Shoulder SP-26 yellow grey           | 15 g    |
| 620127      | ZI-CT/Shoulder SP-27 ivory                 | 15 g    |
| 620128      | ZI-CT/Shoulder SP-28 opaque yellow         | 15 g    |
| 620129      | ZI-CT/Shoulder SP-29 honey                 | 15 g    |
| 620131      | ZI-CT/Opaque Dentine OD-32 havanna         | 15 g    |
| 620132      | ZI-CT/Opaque Dentine OD-37 curry           | 15 g    |
| 620133      | ZI-CT/Opaque Dentine OD-41 orange          | 15 g    |
| 620134      | ZI-CT/Opaque Dentine OD-43 ivory           | 15 g    |
| 620135      | ZI-CT/Opaque Dentine OD-44 cuba            | 15 g    |
| 620146      | ZI-CT/Approximal Dentine AD-1 light yellow | 15 g    |
| 620147      | ZI-CT/Approximal Dentine AD-2 orange       | 15 g    |
| 620151      | ZI-CT/Neck Transpa HT-51 iridescent        | 15 g    |
| 620152      | ZI-CT/Neck Transpa HT-52 kaki              | 15 g    |
| 620153      | ZI-CT/Neck Transpa HT-53 sand              | 15 g    |
| 620154      | ZI-CT/Neck Transpa HT-54 honey             | 15 g    |
| 620155      | ZI-CT/Neck Transpa HT-55 bordeaux          | 15 g    |
| 620156      | ZI-CT/Neck Transpa HT-56 ochre             | 15 g    |
| 620161      | ZI-CT/Make In MI-61 ivory                  | 15 g    |
| 620162      | ZI-CT/Make In MI-62 lemon                  | 15 g    |
| 620163      | ZI-CT/Make In MI-63 honey                  | 15 g    |
| 620164      | ZI-CT/Make In MI-64 flamingo               | 15 g    |
| 620165      | ZI-CT/Make In MI-65 gold                   | 15 g    |
| 620166      | ZI-CT/Make In MI-66 olive                  | 15 g    |
|             | PORCELAIN 20 g                             |         |
| 620001      | ZI-CT/Dentine D-A1                         | 20 g    |
| 620002      | ZI-CT/Dentine D-A2                         | 20 g    |
| 620003      | ZI-CT/Dentine D-A3                         | 20 g    |
| 620004      | ZI-CT/Dentine D-A3,5                       | 20 g    |
| 620005      | ZI-CT/Dentine D-A4                         | 20 g    |
| 620006      | ZI-CT/Dentine D-B1                         | 20 g    |

CREATION ZI-CT | ZIRCONIUM DIOXIDE CERAMIC

| Article no. | Article name                            | Content |
|-------------|-----------------------------------------|---------|
| 620007      | ZI-CT/Dentine D-B2                      | 20 g    |
| 620008      | ZI-CT/Dentine D-B3                      | 20 g    |
| 620009      | ZI-CT/Dentine D-B4                      | 20 g    |
| 620010      | ZI-CT/Dentine D-C1                      | 20 g    |
| 620011      | ZI-CT/Dentine D-C2                      | 20 g    |
| 620012      | ZI-CT/Dentine D-C3                      | 20 g    |
| 620013      | ZI-CT/Dentine D-C4                      | 20 g    |
| 620014      | ZI-CT/Dentine D-D2                      | 20 g    |
| 620015      | ZI-CT/Dentine D-D3                      | 20 g    |
| 620016      | ZI-CT/Dentine D-D4                      | 20 g    |
| 620017      | ZI-CT/Enamel E-57                       | 20 g    |
| 620018      | ZI-CT/Enamel E-58                       | 20 g    |
| 620019      | ZI-CT/Enamel E-59                       | 20 g    |
| 620020      | ZI-CT/Enamel E-60                       | 20 g    |
| 620021      | ZI-CT/Clear CL-O                        | 20 g    |
| 620022      | ZI-CT/Window UC                         | 20 g    |
| 620023      | ZI-CT/Bleach Dentine BD-A               | 20 g    |
| 620024      | ZI-CT/Bleach Dentine BD-B               | 20 g    |
| 620025      | ZI-CT/Bleach Dentine BD-B0              | 20 g    |
| 620026      | ZI-CT/Bleach Enamel S-AB                | 20 g    |
| 620031      | ZI-CT/Effect Enamel SI-01 light yellow  | 20 g    |
| 620032      | ZI-CT/Effect Enamel SI-02 medium yellow | 20 g    |
| 620033      | ZI-CT/Effect Enamel SI-03 heavy yellow  | 20 g    |
| 620034      | ZI-CT/Effect Enamel SI-04 light orange  | 20 g    |
| 620035      | ZI-CT/Effect Enamel SI-05 medium orange | 20 g    |
| 620036      | ZI-CT/Effect Enamel SI-06 heavy orange  | 20 g    |
| 620037      | ZI-CT/Effect Enamel SO-10 blue          | 20 g    |
| 620038      | ZI-CT/Effect Enamel SO-11 orange        | 20 g    |
| 620039      | ZI-CT/Transpa NT neutral                | 20 g    |
| 620040      | ZI-CT/Transpa OT opal                   | 20 g    |
| 620041      | ZI-CT/Transpa TI-1 blue                 | 20 g    |
| 620042      | ZI-CT/Transpa TI-2 white                | 20 g    |
| 620043      | ZI-CT/Transpa TI-3 pink                 | 20 g    |
| 620044      | ZI-CT/Transpa TI-4 yellow               | 20 g    |
| 620045      | ZI-CT/Transpa TI-5 grey                 | 20 g    |
| 620050      | ZI-CT/Pearl Enamel PS-0 white           | 20 g    |
| 620051      | ZI-CT/Pearl Enamel PS-1 blue/white      | 20 g    |
| 620052      | ZI-CT/Pearl Enamel PS-2 flamingo        | 20 g    |
| 620053      | ZI-CT/Pearl Enamel PS-3 orange          | 20 g    |
| 620060      | ZI-CT/Correction Powder KM              | 20 g    |
|             |                                         |         |
|             | <b>PORCELAIN 50 g</b>                   |         |
| 620201      | ZI-CT/Dentine D-A1                      | 50 g    |
| 620202      | ZI-CT/Dentine D-A2                      | 50 g    |
| 620203      | ZI-CT/Dentine D-A3                      | 50 g    |
| 620204      | ZI-CT/Dentine D-A3,5                    | 50 g    |
| 620205      | ZI-CT/Dentine D-A4                      | 50 g    |
| 620206      | ZI-CT/Dentine D-B1                      | 50 g    |

CREATION ZI-CT | ZIRCONIUM DIOXIDE CERAMIC

| Article no. | Article name                               | Content |
|-------------|--------------------------------------------|---------|
| 620207      | ZI-CT/Dentine D-B2                         | 50 g    |
| 620208      | ZI-CT/Dentine D-B3                         | 50 g    |
| 620209      | ZI-CT/Dentine D-B4                         | 50 g    |
| 620210      | ZI-CT/Dentine D-C1                         | 50 g    |
| 620211      | ZI-CT/Dentine D-C2                         | 50 g    |
| 620212      | ZI-CT/Dentine D-C3                         | 50 g    |
| 620213      | ZI-CT/Dentine D-C4                         | 50 g    |
| 620214      | ZI-CT/Dentine D-D2                         | 50 g    |
| 620215      | ZI-CT/Dentine D-D3                         | 50 g    |
| 620216      | ZI-CT/Dentine D-D4                         | 50 g    |
| 620217      | ZI-CT/Enamel E-57                          | 50 g    |
| 620218      | ZI-CT/Enamel E-58                          | 50 g    |
| 620219      | ZI-CT/Enamel E-59                          | 50 g    |
| 620220      | ZI-CT/Enamel E-60                          | 50 g    |
| 620221      | ZI-CT/Clear CL-O                           | 50 g    |
| 620222      | ZI-CT/Window UC                            | 50 g    |
| 620223      | ZI-CT/Bleach Dentine BD-A                  | 50 g    |
| 620224      | ZI-CT/Bleach Dentine BD-B                  | 50 g    |
| 620225      | ZI-CT/Bleach Dentine BD-B0                 | 50 g    |
| 620226      | ZI-CT/Bleach Enamel S-AB                   | 50 g    |
| 620231      | ZI-CT/Effect Enamel SI-01 light yellow     | 50 g    |
| 620232      | ZI-CT/Effect Enamel SI-02 medium yellow    | 50 g    |
| 620233      | ZI-CT/Effect Enamel SI-03 heavy yellow     | 50 g    |
| 620234      | ZI-CT/Effect Enamel SI-04 light orange     | 50 g    |
| 620235      | ZI-CT/Effect Enamel SI-05 medium orange    | 50 g    |
| 620236      | ZI-CT/Effect Enamel SI-06 heavy orange     | 50 g    |
| 620237      | ZI-CT/Effect Enamel SO-10 blue             | 50 g    |
| 620238      | ZI-CT/Effect Enamel SO-11 orange           | 50 g    |
| 620239      | ZI-CT/Transpa NT neutral                   | 50 g    |
| 620240      | ZI-CT/Transpa OT opal                      | 50 g    |
| 620241      | ZI-CT/Transpa TI-1 blue                    | 50 g    |
| 620242      | ZI-CT/Transpa TI-2 white                   | 50 g    |
| 620243      | ZI-CT/Transpa TI-3 pink                    | 50 g    |
| 620244      | ZI-CT/Transpa TI-4 yellow                  | 50 g    |
| 620245      | ZI-CT/Transpa TI-5 grey                    | 50 g    |
| 620250      | ZI-CT/Pearl Enamel PS-0 white              | 50 g    |
| 620251      | ZI-CT/Pearl Enamel PS-1 blue/white         | 50 g    |
| 620252      | ZI-CT/Pearl Enamel PS-2 flamingo           | 50 g    |
| 620253      | ZI-CT/Pearl Enamel PS-3 orange             | 50 g    |
| 620261      | ZI-CT/Opaque Dentine OD-32 havanna         | 50 g    |
| 620262      | ZI-CT/Opaque Dentine OD-37 curry           | 50 g    |
| 620263      | ZI-CT/Opaque Dentine OD-41 orange          | 50 g    |
| 620264      | ZI-CT/Opaque Dentine OD-43 ivory           | 50 g    |
| 620265      | ZI-CT/Opaque Dentine OD-44 cuba            | 50 g    |
| 620266      | ZI-CT/Approximal Dentine AD-1 light yellow | 50 g    |
| 620267      | ZI-CT/Approximal Dentine AD-2 orange       | 50 g    |
| 620271      | ZI-CT/Neck Transpa HT-51 iridescent        | 50 g    |
| 620272      | ZI-CT/Neck Transpa HT-52 kaki              | 50 g    |

CREATION ZI-CT | ZIRCONIUM DIOXIDE CERAMIC

| Article no. | Article name                      | Content |
|-------------|-----------------------------------|---------|
| 620273      | ZI-CT/Neck Transpa HT-53 sand     | 50 g    |
| 620274      | ZI-CT/Neck Transpa HT-54 honey    | 50 g    |
| 620275      | ZI-CT/Neck Transpa HT-55 bordeaux | 50 g    |
| 620276      | ZI-CT/Neck Transpa HT-56 ochre    | 50 g    |
|             |                                   |         |
|             | <b>PORCELAIN 250 g</b>            |         |
| 620302      | ZI-CT/Dentine D-A2                | 250 g   |
| 620303      | ZI-CT/Dentine D-A3                | 250 g   |
| 620318      | ZI-CT/Enamel E-58                 | 250 g   |
| 620319      | ZI-CT/Enamel E-59                 | 250 g   |
| 620321      | ZI-CT/Clear CL-O                  | 250 g   |
|             |                                   |         |
|             | <b>LIQUIDS</b>                    |         |
| 000012      | Modelling Liquid high-fusing ML-H | 25 ml   |
| 000013      | Modelling Liquid high-fusing ML-H | 50 ml   |
| 000014      | Modelling Liquid high-fusing ML-H | 250 ml  |
| 000015      | Modelling Liquid high-fusing ML-H | 1000 ml |
| 000010      | Universal & Glaze Liquid UF/GL    | 25 ml   |
| 000011      | Universal & Glaze Liquid UF/GL    | 50 ml   |
| 000025      | Shoulder Liquid SPL               | 25 ml   |
| 000026      | Shoulder Liquid SPL               | 50 ml   |
| 000090      | Moisturing Fluid                  | 25 ml   |

CREATION ZI-F | LOW-FUSING ZIRCONIUM DIOXIDE CERAMIC

| Article no. | Article name                              | Content |
|-------------|-------------------------------------------|---------|
|             | <b>PORCELAIN 4 g – 15 g</b>               |         |
| 601087      | ZI-F/Frame Shade NT                       | 4 g     |
| 601074      | ZI-F/Glaze GL                             | 10 g    |
| 601081      | ZI-F/Frame Shade FS-1                     | 10 g    |
| 601082      | ZI-F/Frame Shade FS-2                     | 10 g    |
| 601083      | ZI-F/Frame Shade FS-3                     | 10 g    |
| 601084      | ZI-F/Frame Shade FS-4                     | 10 g    |
| 601085      | ZI-F/Frame Shade FS-IL (Illusion)         | 10 g    |
| 601086      | ZI-F/Frame Shade FS-G (Gingiva)           | 10 g    |
| 602051      | ZI-F/Bleach Frame Shade FS-AB             | 10 g    |
| 601091      | ZI-F/Neck Dentine HD-01                   | 15 g    |
| 601092      | ZI-F/Neck Dentine HD-02                   | 15 g    |
| 601093      | ZI-F/Neck Dentine HD-03                   | 15 g    |
| 601094      | ZI-F/Neck Dentine HD-04                   | 15 g    |
| 602056      | ZI-F/Bleach Shoulder SP-AB                | 15 g    |
| 602121      | ZI-F/Shoulder SP-21 neutral               | 15 g    |
| 602122      | ZI-F/Shoulder SP-22 flamingo              | 15 g    |
| 602123      | ZI-F/Shoulder SP-23 sand                  | 15 g    |
| 602124      | ZI-F/Shoulder SP-24 gold                  | 15 g    |
| 602125      | ZI-F/Shoulder SP-25 red brown             | 15 g    |
| 602128      | ZI-F/Shoulder SP-28 opaque yellow         | 15 g    |
| 602145      | ZI-F/Gingiva G1 purple                    | 15 g    |
| 602146      | ZI-F/Gingiva G2 dark pink                 | 15 g    |
| 602147      | ZI-F/Gingiva G3 light pink                | 15 g    |
| 602148      | ZI-F/Gingiva G4 flamingo                  | 15 g    |
| 602149      | ZI-F/Gingiva G5 rose                      | 15 g    |
| 602150      | ZI-F/Gingiva G6 dark pink opaque          | 15 g    |
| 602151      | ZI-F/Gingiva GN neutral                   | 15 g    |
| 602157      | ZI-F/Gingiva G7 violet                    | 15 g    |
| 602161      | ZI-F/Make In MI-61 ivory                  | 15 g    |
| 602163      | ZI-F/Make In MI-63 honey                  | 15 g    |
| 602165      | ZI-F/Make In MI-65 gold                   | 15 g    |
| 652132      | ZI-F/Opaque Dentine OD-32 havanna         | 15 g    |
| 652137      | ZI-F/Opaque Dentine OD-37 curry           | 15 g    |
| 652141      | ZI-F/Opaque Dentine OD-41 orange          | 15 g    |
| 652143      | ZI-F/Opaque Dentine OD-43 ivory           | 15 g    |
| 652144      | ZI-F/Opaque Dentine OD-44 cuba            | 15 g    |
| 652151      | ZI-F/Neck Transpa HT-51 iridescent        | 15 g    |
| 652152      | ZI-F/Neck Transpa HT-52 kaki              | 15 g    |
| 652153      | ZI-F/Neck Transpa HT-53 sand              | 15 g    |
| 652154      | ZI-F/Neck Transpa HT-54 honey             | 15 g    |
| 652156      | ZI-F/Neck Transpa HT-56 ochre             | 15 g    |
| 652168      | ZI-F/Approximal Dentine AD-1 light yellow | 15 g    |
| 652169      | ZI-F/Approximal Dentine AD-2 orange       | 15 g    |
|             |                                           |         |
|             | <b>PORCELAIN 20 g</b>                     |         |
| 652001      | ZI-F/Dentine D-A1                         | 20 g    |
| 652002      | ZI-F/Dentine D-A2                         | 20 g    |

CREATION ZI-F | LOW-FUSING ZIRCONIUM DIOXIDE CERAMIC

| Article no. | Article name                           | Content |
|-------------|----------------------------------------|---------|
| 652003      | ZI-F/Dentine D-A3                      | 20 g    |
| 652004      | ZI-F/Dentine D-A3,5                    | 20 g    |
| 652005      | ZI-F/Dentine D-A4                      | 20 g    |
| 652006      | ZI-F/Dentine D-B1                      | 20 g    |
| 652007      | ZI-F/Dentine D-B2                      | 20 g    |
| 652008      | ZI-F/Dentine D-B3                      | 20 g    |
| 652009      | ZI-F/Dentine D-B4                      | 20 g    |
| 652010      | ZI-F/Dentine D-C1                      | 20 g    |
| 652011      | ZI-F/Dentine D-C2                      | 20 g    |
| 652012      | ZI-F/Dentine D-C3                      | 20 g    |
| 652013      | ZI-F/Dentine D-C4                      | 20 g    |
| 652014      | ZI-F/Dentine D-D2                      | 20 g    |
| 652015      | ZI-F/Dentine D-D3                      | 20 g    |
| 652016      | ZI-F/Dentine D-D4                      | 20 g    |
| 652052      | ZI-F/Bleach Dentine BD-A               | 20 g    |
| 652053      | ZI-F/Bleach Dentine BD-B               | 20 g    |
| 652054      | ZI-F/Bleach Dentine BD-B0              | 20 g    |
| 652055      | ZI-F/Bleach Enamel S-AB                | 20 g    |
| 652057      | ZI-F/Enamel E-57                       | 20 g    |
| 652058      | ZI-F/Enamel E-58                       | 20 g    |
| 652059      | ZI-F/Enamel E-59                       | 20 g    |
| 652060      | ZI-F/Enamel E-60                       | 20 g    |
| 652066      | ZI-F/Effect Enamel SI-02 medium yellow | 20 g    |
| 652068      | ZI-F/Effect Enamel SI-04 light orange  | 20 g    |
| 652070      | ZI-F/Effect Enamel SI-06 heavy orange  | 20 g    |
| 602072      | ZI-F/Window UC                         | 20 g    |
| 652073      | ZI-F/Clear CL-O                        | 20 g    |
| 652074      | ZI-F/Transpa NT neutral                | 20 g    |
| 652075      | ZI-F/Transpa OT opal                   | 20 g    |
| 652076      | ZI-F/Transpa TI-1 blue                 | 20 g    |
| 652077      | ZI-F/Transpa TI-2 white                | 20 g    |
| 652078      | ZI-F/Transpa TI-3 pink                 | 20 g    |
| 652079      | ZI-F/Transpa TI-4 yellow               | 20 g    |
| 652080      | ZI-F/Transpa TI-5 grey                 | 20 g    |
| 652081      | ZI-F/Effect Enamel PS-0 white          | 20 g    |
| 652083      | ZI-F/Effect Enamel PS-3 orange         | 20 g    |
| 652085      | ZI-F/Effect Enamel SO-10 blue          | 20 g    |
| 652086      | ZI-F/Effect Enamel SO-11 orange        | 20 g    |
| 652090      | ZI-F/Correction Powder KM              | 20 g    |
| 602174      | ZI-F/Glaze GL                          | 20 g    |
|             | <b>PORCELAIN 50 g</b>                  |         |
| 605056      | ZI-F/Bleach Shoulder SP-AB             | 50 g    |
| 605121      | ZI-F/Shoulder SP-21 neutral            | 50 g    |
| 605122      | ZI-F/Shoulder SP-22 flamingo           | 50 g    |
| 605123      | ZI-F/Shoulder SP-23 sand               | 50 g    |
| 605124      | ZI-F/Shoulder SP-24 gold               | 50 g    |
| 605125      | ZI-F/Shoulder SP-25 red brown          | 50 g    |

CREATION ZI-F | LOW-FUSING ZIRCONIUM DIOXIDE CERAMIC

| Article no. | Article name                           | Content |
|-------------|----------------------------------------|---------|
| 605128      | ZI-F/Shoulder SP-28 opaque yellow      | 50 g    |
| 605145      | ZI-F/Gingiva G1 purple                 | 50 g    |
| 605146      | ZI-F/Gingiva G2 dark pink              | 50 g    |
| 605147      | ZI-F/Gingiva G3 light pink             | 50 g    |
| 605148      | ZI-F/Gingiva G4 flamingo               | 50 g    |
| 605149      | ZI-F/Gingiva G5 rose                   | 50 g    |
| 605150      | ZI-F/Gingiva G6 dark pink opaque       | 50 g    |
| 605151      | ZI-F/Gingiva GN neutral                | 50 g    |
| 605157      | ZI-F/Gingiva G7 violet                 | 50 g    |
| 605161      | ZI-F/Make In MI-61 ivory               | 50 g    |
| 605163      | ZI-F/Make In MI-63 honey               | 50 g    |
| 605165      | ZI-F/Make In MI-65 gold                | 50 g    |
| 655001      | ZI-F/Dentine D-A1                      | 50 g    |
| 655002      | ZI-F/Dentine D-A2                      | 50 g    |
| 655003      | ZI-F/Dentine D-A3                      | 50 g    |
| 655004      | ZI-F/Dentine D-A3,5                    | 50 g    |
| 655005      | ZI-F/Dentine D-A4                      | 50 g    |
| 655006      | ZI-F/Dentine D-B1                      | 50 g    |
| 655007      | ZI-F/Dentine D-B2                      | 50 g    |
| 655008      | ZI-F/Dentine D-B3                      | 50 g    |
| 655009      | ZI-F/Dentine D-B4                      | 50 g    |
| 655010      | ZI-F/Dentine D-C1                      | 50 g    |
| 655011      | ZI-F/Dentine D-C2                      | 50 g    |
| 655012      | ZI-F/Dentine D-C3                      | 50 g    |
| 655013      | ZI-F/Dentine D-C4                      | 50 g    |
| 655014      | ZI-F/Dentine D-D2                      | 50 g    |
| 655015      | ZI-F/Dentine D-D3                      | 50 g    |
| 655016      | ZI-F/Dentine D-D4                      | 50 g    |
| 655052      | ZI-F/Bleach Dentine BD-A               | 50 g    |
| 655053      | ZI-F/Bleach Dentine BD-B               | 50 g    |
| 655054      | ZI-F/Bleach Dentine BD-B0              | 50 g    |
| 655055      | ZI-F/Bleach Enamel S-AB                | 50 g    |
| 655057      | ZI-F/Enamel E-57                       | 50 g    |
| 655058      | ZI-F/Enamel E-58                       | 50 g    |
| 655059      | ZI-F/Enamel E-59                       | 50 g    |
| 655060      | ZI-F/Enamel E-60                       | 50 g    |
| 655066      | ZI-F/Effect Enamel SI-02 medium yellow | 50 g    |
| 655068      | ZI-F/Effect Enamel SI-04 light orange  | 50 g    |
| 655070      | ZI-F/Effect Enamel SI-06 heavy orange  | 50 g    |
| 605072      | ZI-F/Window UC                         | 50 g    |
| 655073      | ZI-F/Clear CL-O                        | 50 g    |
| 655074      | ZI-F/Transpa NT neutral                | 50 g    |
| 655075      | ZI-F/Transpa OT opal                   | 50 g    |
| 655076      | ZI-F/Transpa TI-1 blue                 | 50 g    |
| 655077      | ZI-F/Transpa TI-2 white                | 50 g    |
| 655078      | ZI-F/Transpa TI-3 pink                 | 50 g    |
| 655079      | ZI-F/Transpa TI-4 yellow               | 50 g    |
| 655080      | ZI-F/Transpa TI-5 grey                 | 50 g    |

CREATION ZI-F | LOW-FUSING ZIRCONIUM DIOXIDE CERAMIC

| Article no. | Article name                              | Content |
|-------------|-------------------------------------------|---------|
| 655081      | ZI-F/Effect Enamel PS-0 white             | 50 g    |
| 655083      | ZI-F/Effect Enamel PS-3 orange            | 50 g    |
| 655085      | ZI-F/Effect Enamel SO-10 blue             | 50 g    |
| 655086      | ZI-F/Effect Enamel SO-11 orange           | 50 g    |
| 655132      | ZI-F/Opaque Dentine OD-32 havanna         | 50 g    |
| 655137      | ZI-F/Opaque Dentine OD-37 curry           | 50 g    |
| 655141      | ZI-F/Opaque Dentine OD-41 orange          | 50 g    |
| 655143      | ZI-F/Opaque Dentine OD-43 ivory           | 50 g    |
| 655144      | ZI-F/Opaque Dentine OD-44 cuba            | 50 g    |
| 655151      | ZI-F/Neck Transpa HT-51 iridescent        | 50 g    |
| 655152      | ZI-F/Neck Transpa HT-52 kaki              | 50 g    |
| 655153      | ZI-F/Neck Transpa HT-53 sand              | 50 g    |
| 655154      | ZI-F/Neck Transpa HT-54 honey             | 50 g    |
| 655156      | ZI-F/Neck Transpa HT-56 ochre             | 50 g    |
| 655168      | ZI-F/Approximal Dentine AD-1 light yellow | 50 g    |
| 655169      | ZI-F/Approximal Dentine AD-2 orange       | 50 g    |
|             |                                           |         |
|             | LIQUIDS                                   |         |
| 000016      | Modelling Liquid low-fusing ML-L          | 25 ml   |
| 000017      | Modelling Liquid low-fusing ML-L          | 50 ml   |
| 000018      | Modelling Liquid low-fusing ML-L          | 250 ml  |
| 000010      | Universal & Glaze Liquid UF/GL            | 25 ml   |
| 000011      | Universal & Glaze Liquid UF/GL            | 50 ml   |
| 000025      | Shoulder Liquid SPL                       | 25 ml   |
| 000026      | Shoulder Liquid SPL                       | 50 ml   |
| 000090      | Moisturing Fluid                          | 25 ml   |

CREATION LS/LS PRESS | LITHIUM DISILICATE

| Article no. | Article name            | Content |
|-------------|-------------------------|---------|
|             | INGOTS 3 g              |         |
| 830101      | LS/LS Press, MO-0       | 5 x 3 g |
| 830102      | LS/LS Press, MO-1       | 5 x 3 g |
| 830103      | LS/LS Press, MO-2       | 5 x 3 g |
| 830104      | LS/LS Press LT-A        | 5 x 3 g |
| 830105      | LS/LS Press LT-B        | 5 x 3 g |
| 830106      | LS/LS Press LT-C        | 5 x 3 g |
| 830107      | LS/LS Press LT-D        | 5 x 3 g |
| 830108      | LS/LS Press MT-B00      | 5 x 3 g |
| 830109      | LS/LS Press MT-B0       | 5 x 3 g |
| 830110      | LS/LS Press MT-A1       | 5 x 3 g |
| 830111      | LS/LS Press MT-A2       | 5 x 3 g |
| 830112      | LS/LS Press MT-A3       | 5 x 3 g |
| 830113      | LS/LS Press MT-B1       | 5 x 3 g |
| 830114      | LS/LS Press MT-B2       | 5 x 3 g |
| 830117      | LS/LS Press MT-D2       | 5 x 3 g |
| 830122      | LS/LS Press HT-B00      | 5 x 3 g |
| 830123      | LS/LS Press HT-B0       | 5 x 3 g |
| 830118      | LS/LS Press HT-1        | 5 x 3 g |
| 830119      | LS/LS Press HT-2        | 5 x 3 g |
| 830120      | LS/LS Press HT-3        | 5 x 3 g |
| 830121      | LS/LS Press HT-4        | 5 x 3 g |
|             |                         |         |
|             | PORCELAIN 10 g – 20 g   |         |
| 830001      | LS/Dentine D-A1         | 20 g    |
| 830002      | LS/Dentine D-A2         | 20 g    |
| 830003      | LS/Dentine D-A3         | 20 g    |
| 830004      | LS/Dentine D-A3,5       | 20 g    |
| 830005      | LS/Dentine D-A4         | 20 g    |
| 830006      | LS/Dentine D-B1         | 20 g    |
| 830007      | LS/Dentine D-B2         | 20 g    |
| 830008      | LS/Dentine D-B3         | 20 g    |
| 830009      | LS/Dentine D-B4         | 20 g    |
| 830010      | LS/Dentine D-C1         | 20 g    |
| 830011      | LS/Dentine D-C2         | 20 g    |
| 830012      | LS/Dentine D-C3         | 20 g    |
| 830013      | LS/Dentine D-C4         | 20 g    |
| 830014      | LS/Dentine D-D2         | 20 g    |
| 830015      | LS/Dentine D-D3         | 20 g    |
| 830016      | LS/Dentine D-D4         | 20 g    |
| 830022      | LS/Bleach Dentine BD-A  | 20 g    |
| 830023      | LS/Bleach Dentine BD-B  | 20 g    |
| 830024      | LS/Bleach Dentine BD-B0 | 20 g    |
| 830040      | LS/Opaque Dentine OD-32 | 15 g    |
| 830041      | LS/Opaque Dentine OD-37 | 15 g    |
| 830042      | LS/Opaque Dentine OD-41 | 15 g    |
| 830043      | LS/Opaque Dentine OD-43 | 15 g    |
| 830017      | LS/Enamel E-57          | 20 g    |

CREATION LS/LS PRESS | LITHIUM DISILICATE

| Article no. | Article name                     | Content |
|-------------|----------------------------------|---------|
| 830018      | LS/Enamel E-58                   | 20 g    |
| 830019      | LS/Enamel E-59                   | 20 g    |
| 830020      | LS/Enamel E-60                   | 20 g    |
| 830025      | LS/Enamel Bleach S-AB            | 20 g    |
| 830021      | LS/Clear CL-O                    | 20 g    |
| 830026      | LS/Transpa NT                    | 20 g    |
| 830027      | LS/Transpa OT                    | 20 g    |
| 830028      | LS/Transpa TI-01                 | 20 g    |
| 830029      | LS/Transpa TI-02                 | 20 g    |
| 830030      | LS/Transpa TI-03                 | 20 g    |
| 830031      | LS/Transpa TI-04                 | 20 g    |
| 830032      | LS/Transpa TI-05                 | 20 g    |
| 830033      | LS/Pearl Enamel PS-0             | 20 g    |
| 830035      | LS/Effect Enamel SI-02           | 20 g    |
| 830036      | LS/Effect Enamel SI-04           | 20 g    |
| 830038      | LS/Effect Enamel SO-10           | 20 g    |
| 830045      | LS/Neck Transpa HT-51            | 15 g    |
| 830046      | LS/Neck Transpa HT-52            | 15 g    |
| 830047      | LS/Neck Transpa HT-53            | 15 g    |
| 830048      | LS/Neck Transpa HT-54            | 15 g    |
| 830050      | LS/Make In MI-61                 | 15 g    |
| 830051      | LS/Make In MI-63                 | 15 g    |
| 830052      | LS/Correction Powder KM          | 20 g    |
| 830055      | LS/Glaze GL                      | 10 g    |
|             | LIQUIDS                          |         |
| 000016      | Modelling Liquid low-fusing ML-L | 25 ml   |
| 000017      | Modelling Liquid low-fusing ML-L | 50 ml   |
| 000018      | Modelling Liquid low-fusing ML-L | 250 ml  |
| 000010      | Universal & Glaze Liquid UF/GL   | 25 ml   |
| 000011      | Universal & Glaze Liquid UF/GL   | 50 ml   |
| 000090      | Moisturing Fluid                 | 25 ml   |

CREACOLOR | STAINS AND MODIFIERS

| Article no. | Article name                     | Content |
|-------------|----------------------------------|---------|
| 001101      | Make Up Neo MN-F1 white          | 4 g     |
| 001102      | Make Up Neo MN-F2 eggshell white | 4 g     |
| 001103      | Make Up Neo MN-F3 lemon yellow   | 4 g     |
| 001104      | Make Up Neo MN-F4 mandarine      | 4 g     |
| 001105      | Make Up Neo MN-F5 flamingo       | 4 g     |
| 001106      | Make Up Neo MN-F6 apricot        | 4 g     |
| 001107      | Make Up Neo MN-F7 light brown    | 4 g     |
| 001108      | Make Up Neo MN-F8 medium brown   | 4 g     |
| 001109      | Make Up Neo MN-F9 olive          | 4 g     |
| 001110      | Make Up Neo MN-F10 red           | 4 g     |
| 001111      | Make Up Neo MN-F11 blue          | 4 g     |
| 001112      | Make Up Neo MN-F12 grey          | 4 g     |
| 001113      | Make Up Neo MN-F13 deep blue     | 4 g     |
| 001114      | Make Up Neo MN-FA stain A        | 4 g     |
| 001115      | Make Up Neo MN-FB stain B        | 4 g     |
| 001116      | Make Up Neo MN-FC stain C        | 4 g     |
| 001117      | Make Up Neo MN-FD stain D        | 4 g     |
| 001121      | Make Up Neo Glaze GL-F           | 4 g     |
|             |                                  |         |
| 001151      | In Nova Neo INN-1 melone         | 3 g     |
| 001152      | In Nova Neo INN-2 rose           | 3 g     |
| 001153      | In Nova Neo INN-3 yellow         | 3 g     |
| 001154      | In Nova Neo INN-4 red brown      | 3 g     |
| 001155      | In Nova Neo INN-5 brown          | 3 g     |
| 001156      | In Nova Neo INN-6 olive          | 3 g     |
| 001157      | In Nova Neo INN-7 grey           | 3 g     |
| 001158      | In Nova Neo INN-8 blue           | 3 g     |
| 001159      | In Nova Neo INN-9 white          | 3 g     |
| 001160      | In Nova Neo INN-10 mais          | 3 g     |
| 001161      | In Nova Neo INN-IL Illusion      | 3 g     |
| 001162      | In Nova Neo INN-CRL Crackliner   | 3 g     |
|             |                                  |         |
| 001125      | Make Up Neo Fluid                | 25 ml   |
| 001126      | Make Up Neo Reconditioner Fluid  | 8 ml    |
| 001163      | In Nova Neo Fluid                | 25 ml   |

CREATION LIQUIDS

| Article no. | Article name                      | Content |
|-------------|-----------------------------------|---------|
| 000012      | Modelling Liquid high-fusing ML-H | 25 ml   |
| 000013      | Modelling Liquid high-fusing ML-H | 50 ml   |
| 000014      | Modelling Liquid high-fusing ML-H | 250 ml  |
| 000015      | Modelling Liquid high-fusing ML-H | 1000 ml |
| 000016      | Modelling Liquid low-fusing ML-L  | 25 ml   |
| 000017      | Modelling Liquid low-fusing ML-L  | 50 ml   |
| 000018      | Modelling Liquid low-fusing ML-L  | 250 ml  |
| 000019      | Opaque Liquid high-fusing OL-H    | 25 ml   |
| 000020      | Opaque Liquid high-fusing OL-H    | 50 ml   |
| 000021      | Opaque Liquid high-fusing OL-H    | 250 ml  |
| 000010      | Universal & Glaze Liquid UF/GL    | 25 ml   |
| 000011      | Universal & Glaze Liquid UF/GL    | 50 ml   |
| 000025      | Shoulder Liquid SPL               | 25 ml   |
| 000026      | Shoulder Liquid SPL               | 50 ml   |
| 000090      | Moisturing Fluid                  | 25 ml   |
| 001125      | Make Up Neo Fluid                 | 25 ml   |
| 001126      | Make Up Neo Reconditioner Fluid   | 8 ml    |
| 001163      | In Nova Neo Fluid                 | 25 ml   |

CREATION CERAMIC ACCESSORIES

| Article no. | Article name                            | Content |
|-------------|-----------------------------------------|---------|
| 000131      | Brush Nr. I                             | 1 pc.   |
| 000132      | Brush Nr. II                            | 1 pc.   |
| 000133      | Brush Nr. III                           | 1 pc.   |
| 000134      | Brush Nr. IV                            | 1 pc.   |
| 000135      | Brush Nr. V                             | 1 pc.   |
| 000136      | Brush Kit                               | 1 pc.   |
| 000321      | Brush Creapast                          | 1 pc.   |
| 000060      | Sketch Pad                              | 1 pc.   |
| 839912      | LS/Manual DE-EN                         | 1 pc.   |
| 910360      | ZI-CT/Manual DE-EN                      | 1 pc.   |
| 911020      | CC/Manual DE-EN                         | 1 pc.   |
| 911506      | Ceramic meets Gingiva                   | 1 pc.   |
| 911510      | VC Brochure                             | 1 pc.   |
| 911517      | VC Workbook                             | 1 pc.   |
| 000059      | CreaPen Separating Pen                  | 1 pc.   |
| 000081      | Creation Flüssige firing pad 10 g       | 1 pc.   |
| 000083      | Creation Flüssige firing pad (3 x 10 g) | 1 pack. |

CREATION SHADE GUIDES

| Article no. | Article name                                                                     | Content |
|-------------|----------------------------------------------------------------------------------|---------|
|             | KIT-BASED SHADE GUIDES                                                           |         |
| 000064      | CC/Shade Guide Kit 1 804                                                         | 1 pc.   |
| 000068      | CC/Shade Guide OM 11-17+G                                                        | 1 pc.   |
| 000073      | CC/Shade Guide Kit 2 1508                                                        | 1 pc.   |
| 000087      | CC/Shade Guide Transition Dentine                                                | 1 pc.   |
| 000092      | ZI-CT/Shade Guide Kit 1                                                          | 1 pc.   |
| 000093      | ZI-CT/Shade Guide Kit 2                                                          | 1 pc.   |
| 000076      | ZI-F.LF.TI/Shade Guide Kit 1 (28 Shades)                                         | 1 pc.   |
| 000077      | ZI-F.LF.TI/Shade Guide Kit 2 (26 Shades)                                         | 1 pc.   |
| 000058      | Creation shade tab                                                               | 1 pc.   |
| 000091      | Shade Guide Wheel Empty                                                          | 1 pc.   |
|             | SHADE GUIDES                                                                     |         |
| 000040      | Shade Guide Dentine A-D (16 tabs)                                                | 1 pc.   |
| 000041      | Shade Guide Opaque Dentine A-D (16 tabs)                                         | 1 pc.   |
| 000042      | Shade Guide Transpa CL, UC, NT, OT (4 tabs)                                      | 1 pc.   |
| 000043      | Shade Guide Opaque Dentine Modifier + AD OD-32, 37, 41, 43, 44, AD-1, 2 (7 tabs) | 1 pc.   |
| 000044      | Shade Guide Enamel E57-E60, SI01-SI06, SO10, SO11, PS0-PS3 (16 tabs)             | 1 pc.   |
| 000045      | Shade Guide Transpa Modifier TI1-TI5, HT51-HT56 (11 tabs)                        | 1 pc.   |
| 000046      | Shade Guide Make In MI61-MI66 (6 tabs)                                           | 1 pc.   |
| 000047      | Shade Guide Shoulder SP21-SP29 (9 tabs)                                          | 1 pc.   |
| 000048      | Shade Guide Gingiva G1-G7, GN (8 tabs)                                           | 1 pc.   |
| 000049      | Shade Guide Bleach BDA, BDB, BDB0, EAB, SP-AB (5 tabs)                           | 1 pc.   |

Shades A1–D4 are referring to the Vita® Shade Guide.  
Vita® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Bad Säckingen, Germany.

AQUALINE

| Article no. | Article name                                                                     | Content |
|-------------|----------------------------------------------------------------------------------|---------|
| AMA-001     | Aqualine master, mixing tray for ceramic, anthracite                             | 1 pc.   |
| AMA-002     | Aqualine master, mixing tray for ceramic, white                                  | 1 pc.   |
| AS-001      | Aqualine smart, mixing tray for ceramic, anthracite                              | 1 pc.   |
| AS-002      | Aqualine smart, mixing tray for ceramic, white                                   | 1 pc.   |
| AM-001      | Aqualine mini, mixing tray für Keramik, anthracite                               | 1 pc.   |
| AM-002      | Aqualine mini, tray box with stain tray, with 15 cavities                        | 1 pc.   |
| AS-014      | Aqualine Pool Set                                                                | 1 pc.   |
| AP-010      | Aqualine master, glass plate, light grey                                         | 1 pc.   |
| AP-011      | Aqualine master, glass plate, dark grey                                          | 1 pc.   |
| AS-010      | Aqualine smart, glass plate, light grey                                          | 1 pc.   |
| AS-011      | Aqualine smart, glass plate, dark grey                                           | 1 pc.   |
| AP-013      | Aqualine master, cellulose strips (40 pcs.)                                      | 1 pack. |
| AS-012      | Aqualine smart, cellulose strips (30 pcs.)                                       | 1 pack. |
| AM-010      | Aqualine mini, cellulose strips (30 pcs.)                                        | 1 pack. |
| AMA-018     | Aqualine master, labels, self-adhesive for stain (5 pcs.)                        | 1 pack. |
| AMA-019     | Aqualine master, labels, self-adhesive for ceramic (20 pcs.)                     | 1 pack. |
| AS-016      | Aqualine master/smart, stain kit, made of porcelain (1 stain tray/1 mixing tray) | 1 pc.   |
| AS-015      | Aqualine Pool Set, pool set sponge (5 pcs.)                                      | 1 pack. |
| AC-001      | Aqualine Glass Care (2 x 30 ml)                                                  | 1 set   |
| AC-002      | Aqualine Two in One Build up Liquid (500 ml)                                     | 1 pc.   |

CREATION VC | LIGHT-CURING COMPOSITE

| Article no. | Article name                          | Content |
|-------------|---------------------------------------|---------|
|             | <b>PORCELAIN 2 ml – 5 ml</b>          |         |
| 820005      | VC/Base Opaque O-BO                   | 2.4 ml  |
| 820006      | VC/Cervical Opaque O-CO               | 2.4 ml  |
| 820011      | VC/Opaque O-A1                        | 2.4 ml  |
| 820012      | VC/Opaque O-A2                        | 2.4 ml  |
| 820013      | VC/Opaque O-A3                        | 2.4 ml  |
| 820014      | VC/Opaque O-A3.5                      | 2.4 ml  |
| 820015      | VC/Opaque O-A4                        | 2.4 ml  |
| 820016      | VC/Opaque O-B1                        | 2.4 ml  |
| 820017      | VC/Opaque O-B2                        | 2.4 ml  |
| 820018      | VC/Opaque O-B3                        | 2.4 ml  |
| 820019      | VC/Opaque O-B4                        | 2.4 ml  |
| 820020      | VC/Opaque O-C1                        | 2.4 ml  |
| 820021      | VC/Opaque O-C2                        | 2.4 ml  |
| 820022      | VC/Opaque O-C3                        | 2.4 ml  |
| 820023      | VC/Opaque O-C4                        | 2.4 ml  |
| 820024      | VC/Opaque O-D2                        | 2.4 ml  |
| 820025      | VC/Opaque O-D3                        | 2.4 ml  |
| 820026      | VC/Opaque O-D4                        | 2.4 ml  |
| 820031      | VC/Opaque Dentine OD-32 havanna       | 2.9 ml  |
| 820032      | VC/Opaque Dentine OD-37 curry         | 2.9 ml  |
| 820033      | VC/Opaque Dentine OD-41 orange        | 2.9 ml  |
| 820034      | VC/Opaque Dentine OD-43 ivory         | 2.9 ml  |
| 820035      | VC/Opaque Dentine OD-44 cuba          | 2.9 ml  |
| 820041      | VC/Dentine D-A1                       | 2.9 ml  |
| 820042      | VC/Dentine D-A2                       | 2.9 ml  |
| 820043      | VC/Dentine D-A3                       | 2.9 ml  |
| 820044      | VC/Dentine D-A3.5                     | 2.9 ml  |
| 820045      | VC/Dentine D-A4                       | 2.9 ml  |
| 820046      | VC/Dentine D-B1                       | 2.9 ml  |
| 820047      | VC/Dentine D-B2                       | 2.9 ml  |
| 820048      | VC/Dentine D-B3                       | 2.9 ml  |
| 820049      | VC/Dentine D-B4                       | 2.9 ml  |
| 820050      | VC/Dentine D-C1                       | 2.9 ml  |
| 820051      | VC/Dentine D-C2                       | 2.9 ml  |
| 820052      | VC/Dentine D-C3                       | 2.9 ml  |
| 820053      | VC/Dentine D-C4                       | 2.9 ml  |
| 820054      | VC/Dentine D-D2                       | 2.9 ml  |
| 820055      | VC/Dentine D-D3                       | 2.9 ml  |
| 820056      | VC/Dentine D-D4                       | 2.9 ml  |
| 820057      | VC/Enamel E-57                        | 2.9 ml  |
| 820058      | VC/Enamel E-58                        | 2.9 ml  |
| 820059      | VC/Enamel E-59                        | 2.9 ml  |
| 820060      | VC/Enamel E-60                        | 2.9 ml  |
| 820061      | VC/Pearl Enamel PS-1 blue             | 2.9 ml  |
| 820062      | VC/Pearl Enamel PS-2 flamingo         | 2.9 ml  |
| 820063      | VC / Effect Enamel SI-01 light yellow | 2.9 ml  |
| 820064      | VC/Effect Enamel SI-02 medium yellow  | 2.9 ml  |

CREATION VC | LIGHT-CURING COMPOSITE

| Article no. | Article name                              | Content |
|-------------|-------------------------------------------|---------|
| 820065      | VC/Effect Enamel SI-03 heavy yellow       | 2.9 ml  |
| 820066      | VC/Transpa TI-01 blue                     | 2.9 ml  |
| 820067      | VC/Transpa TI-02 white                    | 2.9 ml  |
| 820068      | VC/Transpa TI-04 yellow                   | 2.9 ml  |
| 820069      | VC/Transpa TI-05 grey                     | 2.9 ml  |
| 820070      | VC/Neck Transpa HT-52 khaki               | 2.9 ml  |
| 820071      | VC/Neck Transpa HT-53 sand                | 2.9 ml  |
| 820072      | VC/Neck Transpa HT-54 honey               | 2.9 ml  |
| 820080      | VC/Flow Dentine BW                        | 2.0 ml  |
| 820081      | VC/Flow Dentine D-A1                      | 2.0 ml  |
| 820082      | VC/Flow Dentine D-A2                      | 2.0 ml  |
| 820083      | VC/Flow Dentine D-A3                      | 2.0 ml  |
| 820084      | VC/Flow Dentine D-A3.5                    | 2.0 ml  |
| 820085      | VC/Flow Dentine D-B1                      | 2.0 ml  |
| 820086      | VC/Flow Dentine D-B2                      | 2.0 ml  |
| 820087      | VC/Flow Dentine D-B3                      | 2.0 ml  |
| 820088      | VC/Flow Dentine HT                        | 2.0 ml  |
| 820089      | VC/Flow Opaque Dentine OD-A2              | 2.0 ml  |
| 820090      | VC/Flow Opaque Dentine OD-A3              | 2.0 ml  |
| 820091      | VC/Flow Enamel EN-L light                 | 2.0 ml  |
| 820092      | VC/Flow Enamel EN-T transparent           | 2.0 ml  |
| 820100      | VC/Gingiva Opaque O-G3 light pink         | 2.4 ml  |
| 820101      | VC/Gingiva Opaque O-G4 flamingo           | 2.4 ml  |
| 820102      | VC/Gingiva Opaque O-G5 rose               | 2.4 ml  |
| 820103      | VC/Gingiva Opaque O-G6 dark pink          | 2.4 ml  |
| 820104      | VC/Gingiva G3 light pink                  | 2.9 ml  |
| 820105      | VC/Gingiva G4 flamingo                    | 2.9 ml  |
| 820106      | VC/Gingiva G5 rose                        | 2.9 ml  |
| 820107      | VC/Gingiva G6 dark pink                   | 2.9 ml  |
| 820110      | VC/Gingiva Modifier GM-1 purple           | 2.4 ml  |
| 820111      | VC/Gingiva Modifier GM-2 dark pink        | 2.4 ml  |
| 820112      | VC/Gingiva Modifier GM-3 light pink       | 2.4 ml  |
| 820113      | VC/Gingiva Modifier GM-5 rose             | 2.4 ml  |
| 820114      | VC/Gingiva Modifier GM-6 dark pink opaque | 2.4 ml  |
| 820115      | VC/Gingiva Modifier GM-7 violet           | 2.4 ml  |
| 820116      | VC/Gingiva Translucent GT light pink      | 2.4 ml  |
| 820120      | VC/MAKE UP NEO MUN-A stain A              | 2.6 ml  |
| 820121      | VC/MAKE UP NEO MUN-B stain B              | 2.6 ml  |
| 820122      | VC/MAKE UP NEO MUN-C stain C              | 2.6 ml  |
| 820123      | VC/MAKE UP NEO MUN-1 white                | 2.6 ml  |
| 820124      | VC/MAKE UP NEO MUN-2 eggshell             | 2.6 ml  |
| 820125      | VC/MAKE UP NEO MUN-3 lemon yellow         | 2.6 ml  |
| 820126      | VC/MAKE UP NEO MUN-4 mandarine            | 2.6 ml  |
| 820127      | VC/MAKE UP NEO MUN-7 light brown          | 2.6 ml  |
| 820128      | VC/MAKE UP NEO MUN-8 medium brown         | 2.6 ml  |
| 820129      | VC/MAKE UP NEO MUN-9 olive                | 2.6 ml  |
| 820130      | VC/MAKE UP NEO MUN-10 red                 | 2.6 ml  |
| 820131      | VC/MAKE UP NEO MUN-11 blue                | 2.6 ml  |

CREATION VC | LIGHT-CURING COMPOSITE

| Article no. | Article name                     | Content         |
|-------------|----------------------------------|-----------------|
| 820132      | VC/MAKE UP NEO MUN-12 grey       | 2.6 ml          |
| 820133      | VC/MAKE UP NEO MUN-IL Illusion   | 2.6 ml          |
| 820134      | VC/MAKE UP NEO MUN-G Gingiva     | 2.6 ml          |
| 820135      | VC/MAKE UP NEO MUN-CL clear      | 5.0 ml          |
| 820136      | VC/MAKE UP NEO MUN-CL clear h.c. | 5.0 ml          |
| 820027      | VC/Translucent h.v.              | 2.4 ml          |
| 820028      | VC/Translucent l.v.              | 2.4 ml          |
|             | ACCESSORIES                      |                 |
| 820150      | VC/Syringe Plunger               | 1 pc.           |
| 820151      | VC/Separator                     | 5.0 ml          |
| 820154      | VC/Metal Primer                  | 5.0 ml          |
| 820155      | VC/Composite Primer              | 3.0 ml          |
| 820156      | VC/Ceramic Primer                | 3.0 ml          |
| 820157      | VC/Disposable Palette            | 2 x 5 pcs.      |
| 820158      | VC/Brush Round, Size 1           | 2 x 5 pcs.      |
| 820159      | VC/Brush Flat Size 1             | 2 x 5 pcs.      |
| 820160      | VC/Brush Nr. 7                   | 1 pc.           |
| 820161      | VC/Brush Holder                  | 1 pc.           |
| 820152      | VC/Diamond Polisher              | 8 g             |
| 820153      | VC/Air Protect                   | 10.0 ml         |
| 820162      | VC/Dispensing tip                | 20 tips + 1 cap |
| 820163      | VC/Mixing Pad                    | 1 pc.           |

Creation Denture Teeth

CREOPAL | CREAPEARL 2



| CREOPAL DENTURE TEETH SETS |                                   |             |       |         |
|----------------------------|-----------------------------------|-------------|-------|---------|
| Article no.                | Article name                      | Tooth shade | Shape | Content |
| 1100305005                 | Creopal upper jaw anterior tooth  | A1          | AU05  | 1 pc.   |
| 1100305007                 | Creopal upper jaw anterior tooth  | A1          | AU07  | 1 pc.   |
| 1100305010                 | Creopal upper jaw anterior tooth  | A1          | AU10  | 1 pc.   |
| 1100305013                 | Creopal upper jaw anterior tooth  | A1          | AU13  | 1 pc.   |
| 1100305016                 | Creopal upper jaw anterior tooth  | A1          | AU16  | 1 pc.   |
| 1100305019                 | Creopal upper jaw anterior tooth  | A1          | AU19  | 1 pc.   |
| 1100305107                 | Creopal lower jaw anterior tooth  | A1          | AL07  | 1 pc.   |
| 1100305110                 | Creopal lower jaw anterior tooth  | A1          | AL10  | 1 pc.   |
| 1100305113                 | Creopal lower jaw anterior tooth  | A1          | AL13  | 1 pc.   |
| 1100305116                 | Creopal lower jaw anterior tooth  | A1          | AL16  | 1 pc.   |
| 1100305207                 | Creopal upper jaw posterior tooth | A1          | PU07  | 1 pc.   |
| 1100305210                 | Creopal upper jaw posterior tooth | A1          | PU10  | 1 pc.   |
| 1100305213                 | Creopal upper jaw posterior tooth | A1          | PU13  | 1 pc.   |
| 1100305216                 | Creopal upper jaw posterior tooth | A1          | PU16  | 1 pc.   |
| 1100305307                 | Creopal lower jaw posterior tooth | A1          | PL07  | 1 pc.   |
| 1100305310                 | Creopal lower jaw posterior tooth | A1          | PL10  | 1 pc.   |
| 1100305313                 | Creopal lower jaw posterior tooth | A1          | PL13  | 1 pc.   |
| 1100305316                 | Creopal lower jaw posterior tooth | A1          | PL16  | 1 pc.   |
| 1100315005                 | Creopal upper jaw anterior tooth  | A2          | AU05  | 1 pc.   |
| 1100315007                 | Creopal upper jaw anterior tooth  | A2          | AU07  | 1 pc.   |
| 1100315010                 | Creopal upper jaw anterior tooth  | A2          | AU10  | 1 pc.   |
| 1100315013                 | Creopal upper jaw anterior tooth  | A2          | AU13  | 1 pc.   |
| 1100315016                 | Creopal upper jaw anterior tooth  | A2          | AU16  | 1 pc.   |
| 1100315019                 | Creopal upper jaw anterior tooth  | A2          | AU19  | 1 pc.   |
| 1100315107                 | Creopal lower jaw anterior tooth  | A2          | AL07  | 1 pc.   |
| 1100315110                 | Creopal lower jaw anterior tooth  | A2          | AL10  | 1 pc.   |
| 1100315113                 | Creopal lower jaw anterior tooth  | A2          | AL13  | 1 pc.   |
| 1100315116                 | Creopal lower jaw anterior tooth  | A2          | AL16  | 1 pc.   |
| 1100315207                 | Creopal upper jaw posterior tooth | A2          | PU07  | 1 pc.   |
| 1100315210                 | Creopal upper jaw posterior tooth | A2          | PU10  | 1 pc.   |
| 1100315213                 | Creopal upper jaw posterior tooth | A2          | PU13  | 1 pc.   |
| 1100315216                 | Creopal upper jaw posterior tooth | A2          | PU16  | 1 pc.   |
| 1100315307                 | Creopal lower jaw posterior tooth | A2          | PL07  | 1 pc.   |
| 1100315310                 | Creopal lower jaw posterior tooth | A2          | PL10  | 1 pc.   |
| 1100315313                 | Creopal lower jaw posterior tooth | A2          | PL13  | 1 pc.   |
| 1100315316                 | Creopal lower jaw posterior tooth | A2          | PL16  | 1 pc.   |
| 1100325005                 | Creopal upper jaw anterior tooth  | A3          | AU05  | 1 pc.   |
| 1100325007                 | Creopal upper jaw anterior tooth  | A3          | AU07  | 1 pc.   |
| 1100325010                 | Creopal upper jaw anterior tooth  | A3          | AU10  | 1 pc.   |
| 1100325013                 | Creopal upper jaw anterior tooth  | A3          | AU13  | 1 pc.   |
| 1100325016                 | Creopal upper jaw anterior tooth  | A3          | AU16  | 1 pc.   |
| 1100325019                 | Creopal upper jaw anterior tooth  | A3          | AU19  | 1 pc.   |
| 1100325107                 | Creopal lower jaw anterior tooth  | A3          | AL07  | 1 pc.   |
| 1100325110                 | Creopal lower jaw anterior tooth  | A3          | AL10  | 1 pc.   |
| 1100325113                 | Creopal lower jaw anterior tooth  | A3          | AL13  | 1 pc.   |
| 1100325116                 | Creopal lower jaw anterior tooth  | A3          | AL16  | 1 pc.   |
| 1100325207                 | Creopal upper jaw posterior tooth | A3          | PU07  | 1 pc.   |
| 1100325210                 | Creopal upper jaw posterior tooth | A3          | PU10  | 1 pc.   |

8283

8485

C R E O P A L   D E N T U R E   T E E T H   S E T S

| Article no. | Article name                      | Tooth shade | Shape | Content |
|-------------|-----------------------------------|-------------|-------|---------|
| 1100435107  | Creopal lower jaw anterior tooth  | D2          | AL07  | 1 pc.   |
| 1100435110  | Creopal lower jaw anterior tooth  | D2          | AL10  | 1 pc.   |
| 1100435113  | Creopal lower jaw anterior tooth  | D2          | AL13  | 1 pc.   |
| 1100435116  | Creopal lower jaw anterior tooth  | D2          | AL16  | 1 pc.   |
| 1100435207  | Creopal upper jaw posterior tooth | D2          | PU07  | 1 pc.   |
| 1100435210  | Creopal upper jaw posterior tooth | D2          | PU10  | 1 pc.   |
| 1100435213  | Creopal upper jaw posterior tooth | D2          | PU13  | 1 pc.   |
| 1100435216  | Creopal upper jaw posterior tooth | D2          | PU16  | 1 pc.   |
| 1100435307  | Creopal lower jaw posterior tooth | D2          | PL07  | 1 pc.   |
| 1100435310  | Creopal lower jaw posterior tooth | D2          | PL10  | 1 pc.   |
| 1100435313  | Creopal lower jaw posterior tooth | D2          | PL13  | 1 pc.   |
| 1100435316  | Creopal lower jaw posterior tooth | D2          | PL16  | 1 pc.   |
| 1100445005  | Creopal upper jaw anterior tooth  | D3          | AU05  | 1 pc.   |
| 1100445007  | Creopal upper jaw anterior tooth  | D3          | AU07  | 1 pc.   |
| 1100445010  | Creopal upper jaw anterior tooth  | D3          | AU10  | 1 pc.   |
| 1100445013  | Creopal upper jaw anterior tooth  | D3          | AU13  | 1 pc.   |
| 1100445016  | Creopal upper jaw anterior tooth  | D3          | AU16  | 1 pc.   |
| 1100445019  | Creopal upper jaw anterior tooth  | D3          | AU19  | 1 pc.   |
| 1100445107  | Creopal lower jaw anterior tooth  | D3          | AL07  | 1 pc.   |
| 1100445110  | Creopal lower jaw anterior tooth  | D3          | AL10  | 1 pc.   |
| 1100445113  | Creopal lower jaw anterior tooth  | D3          | AL13  | 1 pc.   |
| 1100445116  | Creopal lower jaw anterior tooth  | D3          | AL16  | 1 pc.   |
| 1100445207  | Creopal upper jaw posterior tooth | D3          | PU07  | 1 pc.   |
| 1100445210  | Creopal upper jaw posterior tooth | D3          | PU10  | 1 pc.   |
| 1100445213  | Creopal upper jaw posterior tooth | D3          | PU13  | 1 pc.   |
| 1100445216  | Creopal upper jaw posterior tooth | D3          | PU16  | 1 pc.   |
| 1100445307  | Creopal lower jaw posterior tooth | D3          | PL07  | 1 pc.   |
| 1100445310  | Creopal lower jaw posterior tooth | D3          | PL10  | 1 pc.   |
| 1100445313  | Creopal lower jaw posterior tooth | D3          | PL13  | 1 pc.   |
| 1100445316  | Creopal lower jaw posterior tooth | D3          | PL16  | 1 pc.   |
| 1100455005  | Creopal upper jaw anterior tooth  | D4          | AU05  | 1 pc.   |
| 1100455007  | Creopal upper jaw anterior tooth  | D4          | AU07  | 1 pc.   |
| 1100455010  | Creopal upper jaw anterior tooth  | D4          | AU10  | 1 pc.   |
| 1100455013  | Creopal upper jaw anterior tooth  | D4          | AU13  | 1 pc.   |
| 1100455016  | Creopal upper jaw anterior tooth  | D4          | AU16  | 1 pc.   |
| 1100455019  | Creopal upper jaw anterior tooth  | D4          | AU19  | 1 pc.   |
| 1100455107  | Creopal lower jaw anterior tooth  | D4          | AL07  | 1 pc.   |
| 1100455110  | Creopal lower jaw anterior tooth  | D4          | AL10  | 1 pc.   |
| 1100455113  | Creopal lower jaw anterior tooth  | D4          | AL13  | 1 pc.   |
| 1100455116  | Creopal lower jaw anterior tooth  | D4          | AL16  | 1 pc.   |
| 1100455207  | Creopal upper jaw posterior tooth | D4          | PU07  | 1 pc.   |
| 1100455210  | Creopal upper jaw posterior tooth | D4          | PU10  | 1 pc.   |
| 1100455213  | Creopal upper jaw posterior tooth | D4          | PU13  | 1 pc.   |
| 1100455216  | Creopal upper jaw posterior tooth | D4          | PU16  | 1 pc.   |
| 1100455307  | Creopal lower jaw posterior tooth | D4          | PL07  | 1 pc.   |
| 1100455310  | Creopal lower jaw posterior tooth | D4          | PL10  | 1 pc.   |
| 1100455313  | Creopal lower jaw posterior tooth | D4          | PL13  | 1 pc.   |
| 1100455316  | Creopal lower jaw posterior tooth | D4          | PL16  | 1 pc.   |

C R E O P A L   D E N T U R E   T E E T H   S E T S

| Article no. | Article name                     | Tooth shade | Shape | Content |
|-------------|----------------------------------|-------------|-------|---------|
| 1100465005  | Creopal upper jaw anterior tooth | FW2         | AU05  | 1 pc.   |
| 1100465007  | Creopal upper jaw anterior tooth | FW2         | AU07  | 1 pc.   |
| 1100465010  | Creopal upper jaw anterior tooth | FW2         | AU10  | 1 pc.   |
| 1100465013  | Creopal upper jaw anterior tooth | FW2         | AU13  | 1 pc.   |
| 1100465016  | Creopal upper jaw anterior tooth | FW2         | AU16  | 1 pc.   |
| 1100465019  | Creopal upper jaw anterior tooth | FW2         | AU19  | 1 pc.   |
| 1100465107  | Creopal lower jaw anterior tooth | FW2         | AL07  | 1 pc.   |
| 1100465110  | Creopal lower jaw anterior tooth | FW2         | AL10  | 1 pc.   |
| 1100465113  | Creopal lower jaw anterior tooth | FW2         | AL13  | 1 pc.   |
| 1100465116  | Creopal lower jaw anterior tooth | FW2         | AL16  | 1 pc.   |
| 1100475005  | Creopal upper jaw anterior tooth | FW5         | AU05  | 1 pc.   |
| 1100475007  | Creopal upper jaw anterior tooth | FW5         | AU07  | 1 pc.   |
| 1100475010  | Creopal upper jaw anterior tooth | FW5         | AU10  | 1 pc.   |
| 1100475013  | Creopal upper jaw anterior tooth | FW5         | AU13  | 1 pc.   |
| 1100475016  | Creopal upper jaw anterior tooth | FW5         | AU16  | 1 pc.   |
| 1100475019  | Creopal upper jaw anterior tooth | FW5         | AU19  | 1 pc.   |
| 1100475107  | Creopal lower jaw anterior tooth | FW5         | AL07  | 1 pc.   |
| 1100475110  | Creopal lower jaw anterior tooth | FW5         | AL10  | 1 pc.   |
| 1100475113  | Creopal lower jaw anterior tooth | FW5         | AL13  | 1 pc.   |
| 1100475116  | Creopal lower jaw anterior tooth | FW5         | AL16  | 1 pc.   |



C R E O P A L   S H E L L - S E T S

| Article no. | Article name                            | Tooth shade | Shape | Content |
|-------------|-----------------------------------------|-------------|-------|---------|
| 1200375116  | Creopal Shell lower jaw anterior tooth  | B3          | SAL16 | 1 pc.   |
| 1200375210  | Creopal Shell upper jaw posterior tooth | B3          | SPU10 | 1 pc.   |
| 1200375213  | Creopal Shell upper jaw posterior tooth | B3          | SPU13 | 1 pc.   |
| 1200375216  | Creopal Shell upper jaw posterior tooth | B3          | SPU16 | 1 pc.   |
| 1200375219  | Creopal Shell upper jaw posterior tooth | B3          | SPU19 | 1 pc.   |
| 1200375310  | Creopal Shell lower jaw posterior tooth | B3          | SPL10 | 1 pc.   |
| 1200375313  | Creopal Shell lower jaw posterior tooth | B3          | SPL13 | 1 pc.   |
| 1200375316  | Creopal Shell lower jaw posterior tooth | B3          | SPL16 | 1 pc.   |
| 1200375319  | Creopal Shell lower jaw posterior tooth | B3          | SPL19 | 1 pc.   |
| 1200475010  | Creopal Shell upper jaw anterior tooth  | FW5         | SAU10 | 1 pc.   |
| 1200475013  | Creopal Shell upper jaw anterior tooth  | FW5         | SAU13 | 1 pc.   |
| 1200475016  | Creopal Shell upper jaw anterior tooth  | FW5         | SAU16 | 1 pc.   |
| 1200475019  | Creopal Shell upper jaw anterior tooth  | FW5         | SAU19 | 1 pc.   |
| 1200475110  | Creopal Shell lower jaw anterior tooth  | FW5         | SAL10 | 1 pc.   |
| 1200475113  | Creopal Shell lower jaw anterior tooth  | FW5         | SAL13 | 1 pc.   |
| 1200475116  | Creopal Shell lower jaw anterior tooth  | FW5         | SAL16 | 1 pc.   |

C R E O P A L   D E N T U R E   T E E T H - P R E S E N T A T I O N   B O X

| Article no. | Article name                            | Tooth shade | Shape      | Content |
|-------------|-----------------------------------------|-------------|------------|---------|
| 1100000010  | Creopal Presentation Box, blank         |             |            | 1 pc.   |
| 1100000011  | Creopal Presentation Box, incl. 33 sets | A1          | all shapes | 1 pc.   |
| 1100000012  | Creopal Presentation Box, incl. 33 sets | A2          | all shapes | 1 pc.   |
| 1100000013  | Creopal Presentation Box, incl. 33 sets | A3          | all shapes | 1 pc.   |
| 1100000014  | Creopal Presentation Box, incl. 33 sets | A3.5        | all shapes | 1 pc.   |
| 1100000017  | Creopal Presentation Box, incl. 33 sets | B2          | all shapes | 1 pc.   |
| 1100000018  | Creopal Presentation Box, incl. 33 sets | B3          | all shapes | 1 pc.   |

C R E O P A L   A C C E S S O R I E S

| Article no. | Article name                                     | Tooth shade | Shape | Content |
|-------------|--------------------------------------------------|-------------|-------|---------|
| 1200000001  | Creopal Primer                                   |             |       | 1 pc.   |
| 1100000001  | Creopal Shade Guide                              |             |       | 1 pc.   |
| 1000300030  | Tooth cabinet (blank, without inlays, 5 drawers) |             |       | 1 pc.   |
| PEARL00003  | Assortment inlay for 72 anterior sets            |             |       | 1 pc.   |
| PEARL00004  | Assortment inlay for 48 posterior sets           |             |       | 1 pc.   |

C R E A P E A R L 2   D E N T U R E   T E E T H   S E T S

| Article no. | Article name                                           | Tooth shade | Shape | Content |
|-------------|--------------------------------------------------------|-------------|-------|---------|
| 1000106002  | Creapearl 2 upper jaw anterior tooth Typee Vario       | A1          | V61   | 1 pc.   |
| 1000106003  | Creapearl 2 upper jaw anterior tooth Type Vario        | A1          | V63   | 1 pc.   |
| 1000106004  | Creapearl 2 upper jaw anterior tooth Type Vario        | A1          | V65   | 1 pc.   |
| 1000106005  | Creapearl 2 upper jaw anterior tooth Type Vario        | A1          | V68   | 1 pc.   |
| 1000106006  | Creapearl 2 upper jaw anterior tooth Type Vario        | A1          | V69   | 1 pc.   |
| 1000106007  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A1          | A71   | 1 pc.   |
| 1000106008  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A1          | A73   | 1 pc.   |
| 1000106009  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A1          | A75   | 1 pc.   |
| 1000106010  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A1          | A78   | 1 pc.   |
| 1000106011  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A1          | A79   | 1 pc.   |
| 1000106012  | Creapearl 2 lower jaw anterior tooth Type U            | A1          | U1    | 1 pc.   |
| 1000106013  | Creapearl 2 lower jaw anterior tooth Type U            | A1          | U3    | 1 pc.   |
| 1000106014  | Creapearl 2 lower jaw anterior tooth Type U            | A1          | U5    | 1 pc.   |
| 1000106015  | Creapearl 2 lower jaw anterior tooth Type U            | A1          | U7    | 1 pc.   |
| 1000106016  | Creapearl 2 lower jaw anterior tooth Type U            | A1          | U9    | 1 pc.   |
| 1000106017  | Creapearl 2 upper jaw posterior tooth Type Dynamicline | A1          | DL    | 1 pc.   |
| 1000106018  | Creapearl 2 lower jaw posterior tooth Type Dynamicline | A1          | DL    | 1 pc.   |
| 1000106019  | Creapearl 2 upper jaw posterior tooth Type Dynamicline | A1          | DM    | 1 pc.   |
| 1000106020  | Creapearl 2 lower jaw posterior tooth Type Dynamicline | A1          | DM    | 1 pc.   |
| 1000106021  | Creapearl 2 upper jaw posterior tooth Type Dynamicline | A1          | DS    | 1 pc.   |
| 1000106022  | Creapearl 2 lower jaw posterior tooth Type Dynamicline | A1          | DS    | 1 pc.   |
| 1000106023  | Creapearl 2 upper jaw posterior tooth Type Balanceline | A1          | EM    | 1 pc.   |
| 1000106024  | Creapearl 2 lower jaw posterior tooth Type Balanceline | A1          | EM    | 1 pc.   |
| 1000106025  | Creapearl 2 upper jaw posterior tooth Type Balanceline | A1          | ES    | 1 pc.   |
| 1000106026  | Creapearl 2 lower jaw posterior tooth Type Balanceline | A1          | ES    | 1 pc.   |
| 1000106027  | Creapearl 2 upper jaw posterior tooth Type Balanceline | A1          | FM    | 1 pc.   |
| 1000106028  | Creapearl 2 lower jaw posterior tooth Type Balanceline | A1          | FM    | 1 pc.   |
| 1000106029  | Creapearl 2 upper jaw posterior tooth Type Balanceline | A1          | FS    | 1 pc.   |
| 1000106030  | Creapearl 2 lower jaw posterior tooth Type Balanceline | A1          | FS    | 1 pc.   |
| 1000106031  | Creapearl 2 upper jaw posterior tooth Type Gnathosline | A1          | GM    | 1 pc.   |
| 1000106032  | Creapearl 2 lower jaw posterior tooth Type Gnathosline | A1          | GM    | 1 pc.   |
| 1000106033  | Creapearl 2 upper jaw posterior tooth Type Gnathosline | A1          | GS    | 1 pc.   |
| 1000106034  | Creapearl 2 lower jaw posterior tooth Type Gnathosline | A1          | GS    | 1 pc.   |
| 1000116002  | Creapearl 2 upper jaw anterior tooth Type Vario        | A2          | V61   | 1 pc.   |
| 1000116003  | Creapearl 2 upper jaw anterior tooth Type Vario        | A2          | V63   | 1 pc.   |
| 1000116004  | Creapearl 2 upper jaw anterior tooth Type Vario        | A2          | V65   | 1 pc.   |
| 1000116005  | Creapearl 2 upper jaw anterior tooth Type Vario        | A2          | V68   | 1 pc.   |
| 1000116006  | Creapearl 2 upper jaw anterior tooth Type Vario        | A2          | V69   | 1 pc.   |
| 1000116007  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A2          | A71   | 1 pc.   |
| 1000116008  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A2          | A73   | 1 pc.   |
| 1000116009  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A2          | A75   | 1 pc.   |
| 1000116010  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A2          | A78   | 1 pc.   |
| 1000116011  | Creapearl 2 upper jaw anterior tooth Type Anatomic     | A2          | A79   | 1 pc.   |
| 1000116012  | Creapearl 2 lower jaw anterior tooth Type U            | A2          | U1    | 1 pc.   |
| 1000116013  | Creapearl 2 lower jaw anterior tooth Type U            | A2          | U3    | 1 pc.   |
| 1000116014  | Creapearl 2 lower jaw anterior tooth Type U            | A2          | U5    | 1 pc.   |
| 1000116015  | Creapearl 2 lower jaw anterior tooth Type U            | A2          | U7    | 1 pc.   |
| 1000116016  | Creapearl 2 lower jaw anterior tooth Type U            | A2          | U9    | 1 pc.   |

9293

9495

9697

9899

100101

C R E A P E A R L 2   D E N T U R E   T E E T H - P R E S E N T A T I O N   B O X

| Article no. | Article name                                | Tooth shade | Shape      | Content |
|-------------|---------------------------------------------|-------------|------------|---------|
| CR00000001  | Creopal Presentation Box, blank             |             |            | 1 pc.   |
| 1000300011  | Creapearl 2 Presentation Box, incl. 33 sets | A1          | all shapes | 1 pc.   |
| 1000300012  | Creapearl 2 Presentation Box, incl. 33 sets | A2          | all shapes | 1 pc.   |
| 1000300013  | Creapearl 2 Presentation Box, incl. 33 sets | A3          | all shapes | 1 pc.   |
| 1000300014  | Creapearl 2 Presentation Box, incl. 33 sets | A3,5        | all shapes | 1 pc.   |

C R E A P E A R L 2   A C C E S S O R I E S

| Article no. | Article name                                     | Tooth shade | Shape | Content |
|-------------|--------------------------------------------------|-------------|-------|---------|
| 1000300009  | Creapearl 2 Tooth shade guide                    |             |       | 1 pc.   |
| 1000300030  | Tooth cabinet (blank, without inlays, 5 drawers) |             |       | 1 pc.   |
| PEARL00003  | Assortment inlay for 72 anterior sets            |             |       | 1 pc.   |
| PEARL00004  | Assortment inlay for 48 posterior sets           |             |       | 1 pc.   |

Notice

YOUR LOCAL DEALER



DISTRIBUTOR  
Creation Willi Geller International GmbH | Koblacherstraße 3, 6812 Meiningen, Austria  
Phone +43 (0)5522 76784, Fax +43 (0)5522 76784-59 | info@creation-willigeller.com, www.creation-willigeller.com

MANUFACTURER  
KLEMA Dentalprodukte GmbH | Koblacherstraße 3a, 6812 Meiningen, Austria

WILLI GELLER  
*Creation*  
SIMPLY BRILLIANT